




# Proctoscope 1541

Wing Plast's Proctoscope 1541 is a disposable instrument manufactured in a clear environmentally friendly plastic. The Proctoscope has a lamp holder and is ergonomically and functionally designed for the comfort of the user and the patient. The lamp holder fits Wing Plast's Fibre light adapter 1527.

- Stable construction
- Ease of use
- Optimal illumination in the Proctoscopes distal tip
- A smooth transition between the tip of the obturator and the end of the Proctoscope ensures maximum patient comfort

# Proctoscope 1541



## Material

- Polystyrene crystal, Polyethylene HD
- Melting temperature: 100 °C
- Not autoclave able
- Latex free
- Environmentally friendly


## Weight

- Detail: 40 g
- Section package: 1,1 kg
- Transport packaging: 6,5 kg


## Packaging

- Section package: 25 pcs
- Transport packaging: 150 pcs
- Pallet quantity: 2,400 pcs

## Transport Packaging



## Pallet



30-1541 -001-210104

Kvalitet ISO 9001 Miljö ISO 14001 Medical ISO 13485

**Wing Plast** - Quality Medical Instruments

