



SELEFA® Nonwoven clip cap

REF: 270070, 270071

SELEFA® Nonwoven clip cap is produced from spunbond nonwoven. It has an elastic band around the edge and has a universal size, fitting most users. The cap can be used to cover the hair when working in areas with hygiene demands.

Key Features:

- ▶ Spunbond nonwoven, polypropylene
- ▶ Green or blue
- ▶ Universal size, covers the hair
- ▶ With elastic band
- ▶ For areas with hygienic demands







REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
270070	SELEFA® Nonwoven clip cap, white	Ø 61 cm	100/1000
270071	SELEFA® Nonwoven clip cap, blue	Ø 61 cm	100/1000

Sales Unit: Transport Carton

Technical Summary

Regulations




 Not made with natural rubber latex

Material

Product Part	Raw Material & Characteristic
Material composition	100 % spunbond polypropylene
Elastic band	Synthetic rubber band, latex-free

Warnings and Recommendations

Warnings:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	PE-plastic bag, transparent	
Transport Carton	Carton, ECF or TCF quality	330 x 250 x 270
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
270070	N/A		
		6 438129 355024	6 438129 555028

270071

N/A



OneMed Group Oy
Address: Metsäläntie 20
Postal Address: P.O. Box 10
FI-00321 Helsinki, Finland
Phone +358 20 786 6800
www.onemed.com
Domicile Helsinki
Business ID 2039640-1