

# MEDIKRO

3000 & 1000

## High-Precision Volume Calibration Syringes

- ▶ High Accuracy and Precision
- ▶ Robust Aluminium Design
- ▶ Excellent Stability
- ▶ Available with Universal Adapter Kit



**Product Code:** M9474 (3000ml) • M9477 (1000ml)

Medikro calibration syringes are important tools for laboratory quality management. They are easy to use and meet or exceed the latest ATS/ERS 2005 recommendations for accuracy and precision.

Medikro 3000 and 1000 are made of 6061-T6 aluminium alloy to guarantee robustness and a long product lifetime. Their unique hexagonal shape, with full-length padding on each of the six sides, provides excellent stability. The syringes are designed to maintain integrity in the challenging pulmonary function testing environment.

