

Sharpsafe Trolley


Product Code:

44000000

The Sharpsafe Near Patient Sharps Disposal (NPSD) system is a range of trays and trolleys which allow for sharps to be taken directly to the patient and disposed of at the point of use.

The Sharpsafe Trolley is compatible with a range of Sharpsafe containers and connection adaptors. The trolley comes with a range of four trays to carry equipment and required paraphernalia for injections and other procedures.

Complete with anti-static castors and a single breaking castor for added peace of mind.



Product Details

Product Dimensions (mm):	519 x 415 x 889 (W x D x H)
Product Weight:	5.0 kg
Product Weight (with trays):	6.0 kg
Trolley Material:	Mild bright steel, powder coated in BS00A05 Goose Grey satin
Compatible with Sizes:	9L, 13L, 24L & 30L
Safeclip Compatible:	Yes
Castor Size	50mm
Anti-Static Castors:	Yes
Lockable:	Yes
Number of Trays:	Four
Tray Material:	White PVC Gloss Board
Carton Dimensions (mm):	585 x 458 x 953 (W x D x H)
Carton Weight (when full):	7.4kg
Carton Volume:	0.10m ³

FOR MORE INFORMATION:

✉ sharpsafesales@vernagroup.com

🌐 vernacare.com 📞 +44 (0)1495 235800

Vernacare Limited, Block G, Newbridge Road Industrial Estate, Blackwood, Caerphilly, NP12 2YN, UK
© Copyright Verna Group 2022. All rights reserved. GP0146_A24_01