

0.45L

Product Code:

5170XXXX


Please note the last four digits denote the specific aperture type and lid colour.

Our Sharpsafe containers are specifically designed to enable the safe disposal of clinical sharps and minimise the risk of associated injuries.

Our 0.45L has been purpose designed for use in the community setting and fits neatly into a bag or pocket. The 0.45L is specifically shaped to contain the maximum amount of sharps in the minimum amount of space.

Available with different colour lids for easily identifiable waste segregation and a purpose designed handle to aid with safe use and storage.

Sharpsafe's new 5th Generation containers launched in 2022 and were the world's first full range of sustainable single use sharps containers. Engineered using ultra low carbon footprint recycled material showing our commitment to championing change.



Available Apertures



SUPPORTING
A CIRCULAR
ECONOMY



REDUCING
YOUR CARBON
FOOTPRINT



REDUCING VIRGIN
PLASTIC ENTERING
THE SUPPLY CHAIN



MADE USING LOCALLY
SOURCED LOW
CARBON FOOTPRINT
RECYCLED MATERIALS





Product Certification	
ISO 23907-1:2019	✓
UN3291	✓
BSI Kitemark	✓

Site Certification	
ISO 9001 - Quality Management	✓
ISO 14001 - Environmental Management	✓
ISO 45001 - Occupational Health & Safety	✓

Product Details	
Size:	0.45L
Usable Capacity:	0.36L
Usable Volume:	84%
Weight:	83g
Fill Line Height:	136mm
Autoclavable:	Yes
Autoclavable Conditions:	134°C
Product Dimensions (mm):	170 x 105 x 50 (HxWxD)
Aperture Dimensions (mm):	25 x 10 (LxW)
GTIN:	50562715XXXXX*
Lot/Batch Number:	Not Applicable
Expiry Date:	Not Applicable

Accessory Compatibility			
Wall Bracket	n/a	Pedestal Trolley	n/a
Stand	n/a	Sharpsafe Trolley	n/a
NPSD Tray	n/a	XL Trolley	n/a

Logistic Data	
Carton Quantity:	100
Carton Weight (when empty):	905g
Carton Weight (when full):	9.5kg
Carton Dimensions (mm):	346 x 505 x 445 (HxWxD)
Carton Volume:	0.081m ³
Storage Temperature Conditions:	Ambient
Storage Lower Temp Limit:	None
Storage Upper Temp Limit:	None
Sterile/Non Sterile:	Non Sterile

Features and Benefits			
In-Mould Label	✓	Latex Free	✓
Temporary Closure	✓	Low Carbon Footprint Material	✓
Permanent Closure	✓	Leak Resistance	✓
Moulded Lid	✓	Disconnection Features	✓
Fill Line	✓	Safety Shield	n/a
Protected Access***	n/a	Semi Translucent Lid	n/a
Single Use	✓	Carry Handle	✓
Different Lid Variations**	✓	Puncture Resistant	✓

Pallet Information	
Euro Pallet Carton Quantity:	22
Euro Pallet Dimensions (mm):	1200 x 800
Standard Pallet Carton Quantity:	24
Standard Pallet Dimensions (mm):	1200 x 1000

Manufacturing Information	
Made In:	Wales, UK
By:	Vernacare

* Please note the last five digits are specific to the aperture type
 **Different product codes may apply

FOR MORE INFORMATION:

✉ sharpsafesales@vernagroup.com

🌐 vernacare.com 📞 +44 (0)1495 235800