


Post-Surgical Bra



The Dale Post-Surgical Bra offers the soft support patients need following breast and thoracic surgeries.

Key Features

- ◆ Soft, seamless, breathable fabric stretches to accommodate edema without binding or irritating sensitive tissue.
- ◆ Multi-directional stretch provides uniform coverage over patient's chest and back.
- ◆ Front hook & loop closures facilitate inspection of the surgical site and dressing changes.
- ◆ Detachable straps allow for ease of application for caregivers and patients with limited arm movement.
- ◆ Provides maximum comfort throughout recovery, even while sleeping.
- ◆ Can be used to hold surgical dressings in place without tape.
- ◆ Not made with natural rubber latex.



For Use After a Wide Range of Procedures

- ◆ Oncology
- ◆ Mastectomy
- ◆ Lumpectomy
- ◆ Reconstruction
- ◆ Radiation
- ◆ MammoSite®
- ◆ Biopsy
- ◆ Fine needle aspiration
- ◆ Reduction
- ◆ Breast Lift
- ◆ Thoracic

Product #	Description	Quantity
702	Medium fits B-D 86cm - 92cm (34" - 36")	1/box
703	Large fits B-D 92cm - 97cm (36" - 38")	1/box
704	X-Large fits B-D 97cm - 112cm (38" - 44")	1/box
705	XX-Large fits C-E 117cm - 137cm (46" - 54")	1/box
706	XXX-Large fits E-I 137cm - 152cm (54" - 60")	1/box

HCPCS Code A6457



Dale Medical Products, Inc.

508.695.9316
dalemed.com

Post-Surgical Bra

Features	Clinical Advantages
Soft, Seamless, Multi-Directional Stretch Material	<ul style="list-style-type: none">▪ Soft compression during recovery.▪ Stretches to accommodate edema without binding or irritating sensitive tissue.▪ Provides uniform coverage over patient's chest and back.▪ Stretch knit fabric permits surgical site to breathe and won't cause skin irritation.▪ Enhances recovery by increasing patient confidence and sense of security with mobility.▪ Secures dressings without the need for tape or "ace-type" bandages.▪ Comfortable to wear 24 hours a day.▪ Can be worn post radiation treatment.▪ Machine washable.
Quick-Fastening Hook & Loop Front Closure	<ul style="list-style-type: none">▪ Easy application for patient or caregiver.▪ Facilitates inspection of the surgical site and dressing changes.
Detachable Straps	<ul style="list-style-type: none">▪ Avoids painful arm movements with multiple application options.▪ Easy to apply to a patient in the supine position.▪ Comfortable straps prevent shoulder irritation or "digging".
Convenient Sizing	<ul style="list-style-type: none">▪ Reduces inventory and cost with five sizes available: M, L, XL, XXL, XXXL.▪ Ensures proper sizing for full coverage and gentle support.
Not Made with Natural Rubber Latex	<ul style="list-style-type: none">▪ Eliminates patient and caregiver latex reactions

Learn More

