


PRODUKTINFORMATION-TRAKEALTUB



SILIKONISERAD PVC, ORAL/NASAL

- Försedd med BlueLine® röntgenlinje
- Termosensitiv PVC
- Tydliga markeringar
- 15 mm ISO standard konnektor
- Mjuka avfasade kanter vid tubspets



Benämning	Storlek (mm)	Art nr	Avdfp/st
Trakealtub utan kuff	2,5	100/141/025	10
Trakealtub utan kuff	3,0	100/141/030	10
Trakealtub utan kuff	3,5	100/141/035	10
Trakealtub utan kuff	4,0	100/141/040	10
Trakealtub utan kuff	4,5	100/141/045	10
Trakealtub utan kuff	5,0	100/141/050	10
Trakealtub utan kuff	5,5	100/141/055	10
Trakealtub utan kuff	6,0	100/141/060	10
Trakealtub utan kuff	6,5	100/141/065	10
Trakealtub utan kuff	7,0	100/141/070	10