


Enteral Home Syringes

Enteral ISO SAF Home, reusable enteral syringe with female ENFit connection (ISO 80369-3). Low dose connections are fitted to the 1ml and 3ml syringes.



Enteral ISO SAF Home Order Codes

Size	Low Dose	Codes	NHS SC Code	Packaging	Unit of Sale
1ml	•	ELDISO-1H	FTR2105	1pc/pouch	100pcs./box
3ml	•	ELDISO-3H	FTR2102	1pc/pouch	100pcs./box
5ml		EISO-5H	FTR2101	1pc/pouch	100pcs./box
10ml		EISO-10H	FTR2103	1pc/pouch	100pcs./box
20ml		EISO-20H	FTR2104	1pc/pouch	100pcs./box
60ml		EISO-60H	FTR2106	1pc/pouch	50pcs./box

Product Features

- CE ENTERAL/ORAL
- CE MEASUREMENT
- FULL THREAD
- RUBBER GASKET
- DEHP FREE
- LATEX FREE
- UP TO 84 USES

Which Accessories & Related Products Can I Use?

ENFit BOTTLE ADAPTOR
ISO SAF

See page 34

ENTERAL MEDICINE STRAW
ISO SAF

See page 34

MEDICINE BOTTLE ADAPTOR
ISO SAF

See page 34

ENTERAL SYRINGE CAP
ISO SAF

See page 35