

**ENTERAL**  
ISO SAF

## Enteral Single Use Syringes


Enteral ISO SAF, single use enteral syringes with female ENFit connection (ISO 80369-3). Low dose option available in 1ml & 3ml syringes.



### Enteral ISO SAF Order Codes

Size	Low Dose	Codes	NHS SC Code	Packaging	Unit of Sale
1ml	•	ELDISO-1	FWM2427	1pc/pouch	100pcs./box
3ml	•	ELDISO-3	FWM2428	1pc/pouch	100pcs./box
1ml		EISO-1	FTA350	1pc/pouch	100pcs./box
3ml		EISO-3	FTA351	1pc/pouch	100pcs./box
5ml		EISO-5	FWM2430	1pc/pouch	100pcs./box
10ml		EISO-10	FWM2431	1pc/pouch	100pcs./box
20ml		EISO-20	FWM2426	1pc/pouch	100pcs./box
60ml		EISO-60	FWM2429	1pc/pouch	50pcs./box

### Product Features



### Which Accessories & Related Products Can I Use?

