

ISO Saf ACCESSORIES

ENFit BOTTLE ADAPTOR

ENFit Bottle Adaptor

Available in five sizes to fit the majority of medicine bottles.



Order Codes - Singles

Size	Direct Codes	NHS SC Code	Unit Of Sale
12 – 15mm	BAISO1-S	FDC1366	Single blister 100 / box
14 – 16.5mm	BAISO2-S	FDC1367	Single blister 100 / box
17– 20.5mm	BAISO3-S	FDC1368	Single blister 100 / box
19 - 23mm	BAISO35-S	FDC1369	Single blister 100 / box
24 – 26mm	BAISO4-S	FDC1370	Single blister 100 / box

Order Codes - Bulk Pack

Size	Direct Codes	NHS SC Code	Unit Of Sale
12 – 15mm	BAISO1	FDC1135	100pcs.
14 – 16.5mm	BAISO2	FDC1136	100pcs.
17– 20.5mm	BAISO3	FDC1137	100pcs.
19 - 23mm	BAISO35	FDC1349	100pcs.
24 – 26mm	BAISO4	FDC1138	100pcs.

ENTERAL MEDICINE STRAW

Enteral Medicine Straw

Male ENFit connector
(ISO 80369-3)



Order Codes

Non Sterile

Size	Direct Code	NHS SC Code	Packaging	Unit Of Sale
50mm	MSISO-50	FWM2610	1pc/pouch	100 per box

Dash3 - Sterile

50mm	13051050	FDC1146	1pc/pouch	30 per box
100mm	13051100	FDC1123	1pc/pouch	30 per box
150mm	13051150	FDC1145	1pc/pouch	30 per box

MEDICINE BOTTLE ADAPTOR

Medicine Bottle Adaptor

Stepped Bottle Adaptor, 15-27mm, male ENFit connector (ISO 80369-3).



Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
SBAISO	FDC1134	1pc/pouch	50 per box

ENFit® ACCESSORIES

ENTERAL SYRINGE CAP

ISO  Saf

Enteral Syringe Cap

Male ENFit connector
(ISO 80369-3), sterile.



Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
ECISO	FWD276	1pc/pouch	100 per box
ECISO-10	FWD284	10pc/pouch	100 per box

FILTER NEEDLE

ISO  Saf

Filter Needle

Male ENFit connector (ISO 80369-3),
sterile. 5µm.



Order Codes

Dash3

Size	Direct Code	NHS SC Code	Packaging	Unit Of Sale
18G 38mm 5µm	13061038	FDC1147	1pc/pouch	30pcs./box
18G 50mm 5µm	13061050	FDC1148	1pc/pouch	30pcs./box

ENTERAL SYRINGE CAP

ISO  Saf

ENFit Female Cap

Female ENFit connector
(ISO 80369-3).



Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
FCISO-NS	FWD277	200pc/pouch	200pcs./box
FCISO-10	FWD283	10pc/pouch	150pcs./box

Y-PORT ADAPTOR

ISO  Saf

Y-Port Adaptor

Double male ENFit connector
(ISO 80369-3), converts any
ENFit single port feeding
tube to Y-Port, sterile.



Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
YPA-ISO	FWM2432	1pc/pouch	20 per box