

Disposable Specula and Dispensers

Safe, convenient, and affordable.

- KleenSpec® Disposable Specula are made from non-toxic plastic and are available in a variety of sizes
- KleenScoop™ cerumen removal disposable specula are designed for curettage of the outer ear canal and will reach just beyond the anterior prominence of the ear canal in patients over 10 years old
- The dispenser holds all specula sizes within a series

KleenSpec Diagnostic and KleenScoop Cerumen Removal Disposable Specula Disposable Otoscope Specula

52432-UB	2.75 mm Universal KleenSpec® Paediatric Disposable Ear Specula—1 Box, Qty 850
52432-U	2.75 mm Universal KleenSpec® Paediatric Disposable—1 Case, Qty 8500
52434-UB	4.25 mm Universal KleenSpec Disposable Ear Specula—1 Box, Qty 850
52434-U	4.25 mm Universal KleenSpec Disposable Ear Specula—1 Case, Qty 8500
52450	KleenScoop Specula (case of 10 boxes, 100 pieces per box)

Disposable Otoscope Specula Dispenser

52400	Dispenser for Disposable Ear Specula without Storage, 524 Series (case of 10 Dispensers)
52401	Dispenser for Disposable Ear Specula with Storage, 524 Series (case of 10 Dispensers)
52101	Dispenser for Disposable Ear Specula with Storage, 521 Series (case of 10 Dispensers)
52100	Dispenser for Disposable Ear Specula without Storage, 521 Series (case of 10 Dispensers)

KleenSpec Pneumatic, Operating and Consulting Disposable Specula—521 Series Disposable Otoscope Specula

52133-B	3.0 mm KleenSpec Disposable Ear Specula for Pneumatic, Operating and Consulting Otoscopes (case of 10 bags, 500/bag)
52134	4.0 mm KleenSpec Disposable Ear Specula for Pneumatic, Operating and Consulting Otoscopes (case of 10 bags, 500/bag)
52135-B	5.0 mm KleenSpec Disposable Ear Specula for Pneumatic, Operating and Consulting Otoscopes (case of 10 bags, 500/bag)



52401
Disposable Otoscope
Specula Dispenser for 524
Series Otoscope Specula

2.75 mm

4.25 mm



52100
Disposable Otoscope Specula
Dispenser for 521 Series
Otoscope Specula

3.0 mm

4.0 mm

5.0 mm