


GE Healthcare

TruSignal sensors and cables – monitor compatibility

GE connector


CARESCAPE® modular monitors and S/5, AS/3 and CS/3 modular monitors with E-(P)RESTN or E-PSM(P) modules. CARESCAPE V100 with TruSignal® option, Responder, B20 and B40 monitors.



- TruSignal** Interconnect cable with GE connector
TS-G3 3 m/10 ft
- TruSignal** Integrated Finger sensor with GE connector
TS-F2-GE 2 m/6.6 ft
TS-F4-GE 4 m/13 ft
- TruSignal** Integrated FingerTip sensor with GE connector
TS-SA4-GE 4 m/13 ft
- TruSignal** Integrated PediTip sensor with GE connector
TS-SP3-GE 3 m/10 ft
- TruSignal** Integrated Ear sensors with GE connector
TS-E2-GE 2 m/6.6 ft
TS-E4-GE 4 m/13 ft

TruSat connector


GE Ohmeda TruSat® pulse oximeter.



- TruSignal** Interconnect cable with TruSat connector
TS-M3 3 m/10 ft
- TruSignal** Integrated Finger sensor with TruSat connector
TS-F4-MC 4 m/13 ft
- TruSignal** Integrated FingerTip sensor with TruSat connector
TS-SA3-MC 3 m/10 ft
- TruSignal** Integrated PediTip sensor with TruSat connector
TS-SP3-MC 3 m/10 ft
- TruSignal** Integrated Ear sensor with TruSat connector
TS-E4-MC 4 m/13 ft

Ohmeda connector


Handheld pulse oximeters (GE Ohmeda TuffSat®, 3775), tabletop pulse oximeters (3700 series, 3800, 3900, 3900P), respiratory gas monitors (4700, 5250 RGM), CARESCAPE modular monitors and S/5, AS/3, CS/3 modular monitors with M-OSAT module, CardiCap®/5 with N-XOSAT option.




- TruSignal** Interconnect cable with Ohmeda connector
TS-H3 3 m/10 ft
- TruSignal** Integrated Finger sensor with Ohmeda connector
TS-F1-H 1 m/3.3 ft
TS-F4-H 4 m/13 ft
- TruSignal** Integrated Ear sensor with Ohmeda connector
TS-E4-H 4 m/13 ft
- TruSignal** Integrated Wrap sensor with Ohmeda connector
TS-W4-H 4 m/13 ft

Datex connector

CARESCAPE modular monitors and S/5, AS/3 and CS/3 modular monitors with M-NE12ST(P)R, M-(N)EST(P)(R) or M-(P)RESTN modules; CardiCap®/5 and S/5 Light Monitors.



- TruSignal** Interconnect cable with Datex connector
TS-N3 3 m/10 ft
- TruSignal** Integrated Finger sensor with Datex connector
TS-F4-N 4 m/13 ft
- TruSignal** Integrated Ear sensor with Datex connector
TS-E4-N 4 m/13 ft



Additional information on monitor compatibility is available upon request.
© 2013 General Electric Company – All rights reserved.
GE and GE monogram are trademarks of General Electric Company.
*Trademark of General Electric Company.
General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.






















TruSignal sensors and cables – patient connections

(For more information refer to the sensor’s instructions for use)

Patient weight ranges and application suggestions

Preferably apply the sensor to a site where the light source and detector are aligned

Sensor selection considerations

TruSignal Adult/Pediatric sensor				20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)	Older pediatric through adult application												
	TS-AP-10 TS-AP-25	10/box 25/box				- Medium- to long-term monitoring - Recommended for: - high-motion environments - infection control - low-perfusion conditions												
TruSignal Adult Adhesive Wrap sensor				>20 kg (>44 lb)		Older pediatric through adult application												
	TS-AAW-10 TS-AAW-25	10/box 25/box				- Medium- to long-term monitoring - Recommended for: - high-motion environments - infection control - low-perfusion conditions												
TruSignal Pediatric Adhesive Wrap sensor				3-20 kg (6.6-44 lb)		Infant through older pediatric application												
	TS-PAW-10 TS-PAW-25	10/box 25/box				- Medium- to long-term monitoring - Recommended for: - high-motion environments - infection control - low-perfusion conditions - patients with fragile skin												
TruSignal AllFit Sensor		Attachment Accessories		≤ 3 kg (≤ 6.6 lb)	3-20 kg (6.6-44 lb)	20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)	Infant through adult application										
	TS-AF-10 TS-AF-25	10/box 25/box	Replacement tape included		(TruTrak/TruTrak+* and TruSignal monitors only)					- Medium- to long-term monitoring - Recommended for: - high-motion environments - infection control - patients with long fingernails - arthritic patients								
TruSignal Sensitive Skin Sensor		Attachment Accessories		≤ 3 kg (≤ 6.6 lb)	3-20 kg (6.6-44 lb)	20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)	Premature infant/neonate through small adult application										
	TS-SE-3	3/box	Foam wrap and adhesive tape included		OXY-SND Infant foam sandal	3/box		OXY-RWS, OXY-RWM Replacement small and medium foam wrap	24/box		OXY-RTW Replacement wide adhesive tape	100/box	(TruTrak/TruTrak+ and TruSignal monitors only)					- Medium- to long-term monitoring - Recommended for: - high-motion environments if used with adhesive tape - fragile skin if used with foam wrap or infant foam sandal - patients with long fingernails
TruSignal Wrap Sensor		Attachment Accessories		3-20 kg (6.6-44 lb)		20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)	Infant/neonate through adult application										
	TS-W-D, TS-W4-H	1/box	Foam wrap and adhesive tape included		OXY-RWL Replacement large foam wrap	24/box		OXY-RTW Replacement wide adhesive tape	100/box				- Short- to medium-term monitoring - Recommended for: - high-motion environments if used with adhesive tape - fragile skin if used with foam wrap or infant foam sandal - patients with long fingernails					
TruSignal Finger Sensor				≥ 20 kg (≥ 44 lb)		Older pediatric through adult application												
	TS-F-D, TS-F2-GE, TS-F4-GE, TS-F4-MC, TS-F1-H, TS-F4-H, TS-F4-N	1/box					- Short-term monitoring - Recommended for: - spot-checks											
TruSignal FingerTip Sensor				> 30 kg (> 66 lb)		Older pediatric through adult application												
	TS-SA-D, TS-SA4-GE, TS-SA3-MC	1/box				Short- to long-term monitoring												
TruSignal PediTip Sensor				15-30 kg (33-66 lb)		Pediatric through small adult application												
	TS-SP-D, TS-SP3-GE, TS-SP3-MC	1/box				Short- to long-term monitoring												
TruSignal Ear Sensor		Attachment Accessories		≥ 10 kg (≥ 22 lb)		Pediatric through adult application												
	TS-E-D, TS-E2-GE, TS-E4-GE, TS-E4-MC, TS-E4-N, TS-E4-H	1/box	Headband and bedsheet clip included		OXY-HB Replacement Headband	3/box		- Short- to medium-term monitoring - Recommended for: - use during stress testing with use of headband - low-perfusion conditions - poor peripheral circulation										