

ETHILON*

Monofilament Polyamid 6

SUTUR

Beskrivelse:

Monofilament

Syntetisk

Ikke absorberbar

Trenger lett gjennom vevet.

Meget vennlig mot vevet.

Opprettholder permanent styrke i vevet.

Farge:

Blå, svart eller gjennomsiktig.

Variasjonsbredde:

11-0 til 2 (USP). Leveres som sutur med nål, for eksempel med nåler fra MULTIPASS og PRIME. ETHILON* er også å få som ligatur.

Indikasjon

ETHILON*-sutur (svart eller gjennomsiktig) er indikert til generelt bløtt vev og/eller til underbinding, herunder til bruk i forbindelse med kardiovaskulære, oftamologiske, mikrokirurgiske og nevrokirurgiske inngrep. Blå ETHILON* sutur er indikert til suturering av hud. Typiske bruksområder er som følger:

- Suturering av hud
- Nervereparasjon
- Oftamologi

Nærmere opplysninger i bruksanvisningen som følger med hver pakning.



ETHILON*

Monofilament Polyamid 6

ETHICON
a Johnson & Johnson company


THE QUALITY
YOUR PATIENTS DESERVE


ETHILON*
Monofilament Polyamide 6

3/8 Circle Reverse Cutting PRIME® Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-6	8 mm	7-0	45		36	1647H
 P-1	11 mm	7-0	45		36	1696H
		7-0	45	MULTIPASS*	36	MPE1696H
		6-0	45		36	697H
		6-0	45	MULTIPASS*	36	MPE697H
 P-3	13 mm	6-0	45		36	1698H
		6-0	45	MULTIPASS*	36	MPE1698H
		5-0	45		36	698H
		5-0	45		36	EH7800H
		5-0	45	MULTIPASS*	36	MPE698H
		5-0	45	MULTIPASS*	36	MPE7800H
		5-0	45		24	W1852T
		4-0	45		36	691H
		4-0	45		36	699H
		4-0	45		36	EH7801H
		4-0	45	MULTIPASS*	36	MPE691H
		4-0	45	MULTIPASS*	36	MPE699H
		4-0	45	MULTIPASS*	36	MPE7801H
		4-0	45	MULTIPASS*	36	MPE7862H
 PS-3	16 mm	6-0	45		24	W1614T
		5-0	45		36	EH7930H
		5-0	45	MULTIPASS*	36	MPE7930H
		4-0	45		36	EH7931H
		4-0	45	MULTIPASS*	36	MPE7931H
		4-0	45		24	W1861T
		5-0	45		36	EH7804H
 PS-2	19 mm	5-0	45	MULTIPASS*	36	MPE7804H
		4-0	45		36	EH7162H
		4-0	45		24	W1620T
		4-0	45	MULTIPASS*	36	MPE1667H
		4-0	45	MULTIPASS*	36	MPE1855H
		3-0	45		36	EH7803H
		3-0	45	MULTIPASS*	36	MPE1621H
		4-0	45	SLIMMED	36	1667H
 PS-1	24 mm	4-0	45		36	EH7957H
		4-0	75		36	EH7956H
		4-0	75	MULTIPASS*	36	MPE7956H
		3-0	45		24	W1857T
		3-0	75		36	EH7954H
		3-0	75	MULTIPASS*	36	MPE7954H

ETHILON*

Monofilament Polyamid 6


3/8 Circle Reverse Cutting PRIME [®] Needles						
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box
 PS	26 mm	3-0	●	45	MULTIPASS [®]	36
		3-0	●	45		24
		3-0	●	75		24
 PSL	30 mm	3-0	●	75		36
 PSLX	36 mm	3-0	●	75		24

3/8 Circle Reverse Cutting Needles						
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box
 G-2	13 mm	6-0	●	45		36
		5-0	●	45		36
 M-2	13 mm	5-0	●	45		36
 FS-3	16 mm	6-0	●	35		12
		6-0	●	45	VISI-BLACK [®]	36
		5-0	●	45		36
		5-0	●	45		36
		4-0	●	45		36
 FS-2	19 mm	5-0	●	45		36
		4-0	●	45		36
		4-0	●	75		36
		3-0	●	45		36
 FS-2S	19 mm	4-0	○	45	SLIMMED	36
		4-0	●	45	SLIMMED	36
 FS-1	24 mm	5-0	●	45		36
		4-0	●	45		36
		4-0	●	45		36
		4-0	●	75		36
		3-0	●	45		36
		3-0	●	75		36
		2-0	●	45		36
		2-0	●	75		36
 FS	26 mm	4-0	●	45		36
		3-0	●	45		36
		2-0	●	45		36
 FSL	30 mm	4-0	●	75		36
		3-0	●	45		36
		3-0	●	75		36
		2-0	●	75		36

ETHILON*

Monofilament Polyamid 6

3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FSLX 36 mm	3-0	●	75		36	1673H
	2-0	●	75		36	1674H
	0	●	75		36	EH7809H
 LR-60 60 mm	2-0	●	100		12	W558
 LR 76 mm	2	●	3x50		12	3460TG


3/8 Circle Conventional Cutting PRIME® Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-1 13 mm	6-0	●	45		36	JM5140H
 PC-3 16 mm	6-0	●	45	MULTIPASS®	36	MPE1853H
	6-0	●	45		24	W1615T
	6-0	●	45		24	W1853T
	5-0	●	45	MULTIPASS®	36	MPE1616H
	5-0	●	45		24	W1616T
	5-0	●	45		24	W1854T
 PC-5 19 mm	4-0	●	45	MULTIPASS®	36	MPE1856H
	4-0	●	45		24	W1856T
 PC-25 26 mm	6-0	●	45		24	W1859T
	5-0	●	45		24	W1858T
	5-0	●	45		12	W532
 KP-3 90 mm	2-0	●	100		12	W795
	1	●	100		12	W797
	1	●	2x50	With retention tubing	12	W2797
 KP-3 90 mm	1	●	75		12	W689


ETHILON*

Monofilament Polyamid 6

▼ 1/2 Circle Reverse Cutting Needles (Heavier Wire Diameter)						
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box
 M-1	13 mm	5-0	●	45		36
 J-1	17 mm	5-0	●	75		36
 X-1	22 mm	3-0	●	45		36
 CPX	48 mm	1	●	100		12
						Code
						648H
						F2423
						EH7943H
						W738

▲ Straight Cutting Needles						
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box
 KS	60 mm	3-0	●	75		36
		2-0	●	75		36
						Code
						627H
						628H

◎ 1/2 Circle Taper Point Needles						
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box
 MH-1	31 mm	0	●	150	Loop	12
						Code
						W739

◎ 1/2 Circle Taper Point Needles (Heavier Wire Diameter)						
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box
 CT	40 mm	0	●	200	Loop	12
 CTX	48 mm	1	●	150	Loop	12
						Code
						W760
						W748

◎ Straight Taper Point Needles						
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box
 STAW	65 mm	2-0	●	75		24
						Code
						EH7625E

ETHILON*

Monofilament Polyamid 6

ETHICON
a Johnson & Johnson company

THE QUALITY
YOUR PATIENTS DESERVE

ETHILON*
Monofilament Polyamide 6



3/8 Circle Taper Point Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 BV 140-3	3,7 mm	10-0	●	15	Non-sterile	12	NS2820
 BV130-3	3,7 mm	9-0	●	13		12	FG2819
 BV75-3	3,8 mm	10-0	●	13		12	W2870
		9-0	●	13		12	W2871
 BV50-3	3,8 mm	11-0	●	13		12	W2881
 BV130-4	4,7 mm	9-0	●	13		12	W2813
		8-0	●	13		12	W2812
 BVH100-3	5 mm	10-0	●	13		12	W2860
 BV75-4	5,1 mm	10-0	●	13		12	W2850
 BV-6	5,1 mm	11-0	●	10		12	EH7438G
 BV-4	5,1 mm	10-0	●	15		12	2814G
 BV100-4	5,1 mm	10-0	●	13		12	W2830
		9-0	●	13		12	W2829
		9-0	●	15		12	EH7448G
 BV-2	6,5 mm	10-0	●	15		12	2810G
		9-0	●	15		12	2809G
		8-0	●	15		12	2808G
 BV130-5	6,5 mm	10-0	●	13		12	W2810
		8-0	●	13		12	W2808
 BV130-5	6,5 mm	10-0	●	30		12	W2801
 BV175-6	8 mm	8-0	●	45		12	W8170
 BV-1	9,3 mm	7-0	●	45		36	EH7446H
 C-1	13 mm	6-0	●	75		36	EH7500H



3/8 Circle CS ULTIMA® Spatula

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
UCS-35	6,2 mm	10-0	●	15	12	7060G
CS-140-6	6,5 mm	10-0	●	30	12	U7003
CS-140-6	6,5 mm	10-0	●	30	12	U7061

ETHILON*

Monofilament Polyamid 6

3/8 Circle MICROPOINT® Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-35-C 	6 mm	10-0	●	15	12	EH7980G
GS-9  	6,5 mm	10-0	●	30	12	1718G
		9-0	●	30	12	1717G
TG140-6  	6,5 mm	11-0	●	30	12	W1780
		10-0	●	30	12	W1770
		9-0	●	30	12	W1769

1/2 Circle MICROPOINT® Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
G-7 	8 mm	8-0	●	45	12	W1765

1/2 Circle MICROPOINT® Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG160-6 	5,5 mm	10-0	●	30	12	W1756
TG175-8 	7 mm	8-0	●	30	12	1716G

1/2 Circle CS ULTIMA® Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-160-6 	5,5 mm	10-0	●	30	12	U7000

1/4 Circle MICROPOINT® Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG100-8 	6,5 mm	8-0	●	30	12	B7714

1/4 Circle Spatula Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
S-14 	8 mm	5-0	●	20	36	MN5132H

ETHILON*

Monofilament Polyamid 6



THE QUALITY
YOUR PATIENTS DESERVE

ETHILON*
Monofilament Polyamide 6

 1/4 Circle TAPERCUT® Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V100-3	3,5 mm	10-0	●	13		12 W2899
		9-0	●	13		12 W2898
 V150-4	4,5 mm	10-0	●	15		12 W1701

 Compound Curve CS ULTIMA® Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CS-C-6	5 mm	10-0	●	30		12 U7009

Ligatures						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
SUTUPAK®	4-0	●	3x45		36	EH8751H
	3-0	●	10x45		36	EH8562H
LIGAPAK®	2-0	●	250	Reel	20	F4354