

VICRYL*

(Polyglactin 910) Suture



VICRYL*

(Polyglactin 910) Braided and Monofilament

SUTURE

Description:	Braided	For easy handling and secure knot tying.
	Coated	For smooth passage through tissue and easy knot tie down.
	Synthetic	For minimal tissue reaction.
	Absorbable	By hydrolysis. Essentially complete in 56-70 days.

Tensile strength:	Post implantation	Approximate % original strength remaining
	14 days	75%
	21 days	50% (6-0 and larger)
		40% (7-0 and smaller)
	28 days	25% (6-0 and larger)

Color: Violet or undyed.

Range: 10-0 to 6 (USP). Supplied as needled sutures and ligatures. Sizes 10-0 and 9-0 are monofilament in structure.

Indications: VICRYL* sutures are intended for use in soft tissue approximation and/or ligation, including use in ophthalmic surgery, peripheral nerve adaptation and microsurgery for vessels less than 2 mm in diameter. Typical areas of use include:








- Fascia closure - due to the 28-day claim
- Subcutaneous fat
- Joint capsule
- Uterus

Full details in the Instruction for Use included in every package.







VICRYL* (Polyglactin 910) Braided or Monofilament Suture

1/2 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 TF-6	6,5 mm	7-0	VIOLET	30		36	V1028H
 TF-1	10 mm	5-0	VIOLET	70		36	V133H
 RB-2	13 mm	6-0	VIOLET	70		36	PZ2571H
		6-0	UNDYED	45		12	W9981
		5-0	UNDYED	70	450 microns	36	V433H
 TF	13 mm	6-0	VIOLET	70		36	V926H
		5-0	VIOLET	70	350 microns	36	V925H
		4-0	VIOLET	70		36	V924H
 RB-1	17 mm	6-0	VIOLET	70		36	V302H
		5-0	VIOLET	70		36	V303H
		5-0	VIOLET	70	VISI-BLACK*	36	V3030H
		5-0	UNDYED	70		36	V213H
		4-0	VIOLET	70		36	V304H
		4-0	VIOLET	70	VISI-BLACK*	36	V3040H
		4-0	UNDYED	70		36	V214H
		4-0	VIOLET	8x45 C/R		12	V734D
		4-0	VIOLET	4x45 C/R		12	V735G
		4-0	VIOLET	8x70 C/R		12	V799D
		3-0	VIOLET	70		36	V305H
		3-0	UNDYED	70		36	V215H
		3-0	VIOLET	8x45 C/R		12	J713D
2-0	VIOLET	70		36	V306H		
 RB-1 RB-1	17 mm	3-0	VIOLET	70		36	V667H
 SH-1	22 mm	4-0	VIOLET	70		36	V310H
		4-0	VIOLET	70	VISI-BLACK*	36	V3100H
		4-0	UNDYED	70		36	V426H
		3-0	VIOLET	70		36	V311H
		3-0	VIOLET	70	VISI-BLACK*	36	V3110H
		3-0	UNDYED	70		36	V219H
		3-0	VIOLET	4x45 C/R		36	V782H
		2-0	VIOLET	70		36	V312H
		2-0	VIOLET	70	VISI-BLACK*	36	V3120H
		2-0	UNDYED	70		36	V220H
		2-0	VIOLET	8x45 C/R		24	V792E









VICRYL* (Polyglactin 910) Braided or Monofilament Suture

1/2 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 SH	26 mm	4-0	VIOLET	70		36	V315H
		4-0	VIOLET	70	VISI-BLACK*	36	V3150H
		4-0	VIOLET	8x70 C/R		12	V783G
		3-0	VIOLET	70		36	V316H
		3-0	VIOLET	70	VISI-BLACK*	36	V3160H
		3-0	UNDYED	70		36	V416H
		3-0	UNDYED	70	VISI-BLACK*	36	V4416H
		3-0	VIOLET	90		36	V242H
		3-0	VIOLET	8x45 C/R		12	J774D
		3-0	VIOLET	8x45 C/R	VISI-BLACK*	24	V7740E
		3-0	VIOLET	8X70 C/R	VISI-BLACK*	24	V7840E
		2-0	VIOLET	70		36	V317H
		2-0	VIOLET	70	VISI-BLACK*	36	V3170H
		2-0	UNDYED	70		36	V224H
		2-0	VIOLET	90		36	V243H
		2-0	VIOLET	4X45 C/R		36	V779H
		2-0	VIOLET	8x45 C/R		12	J775D
		0	VIOLET	70		36	V318H
0	UNDYED	70		36	V225H		
 MH-1	31 mm	3-0	VIOLET	70		36	V245H
		2-0	VIOLET	70		36	V246H
		0	VIOLET	70		36	V247H
 MH	36 mm	4-0	VIOLET	90	Non-stocked item	36	PG2592H
		3-0	VIOLET	70		36	V322H
		3-0	VIOLET	90		36	V100H
		2-0	VIOLET	70		36	V323H
		2-0	UNDYED	70		36	V226H
		2-0	VIOLET	90		36	V101H
		0	VIOLET	70		36	V324H
		0	VIOLET	90		36	V102H
 LH	40 mm	2-0	VIOLET	75		12	W9150






VICRYL* (Polyglactin 910) Braided or Monofilament Suture

1/2 Circle Taper Point Needles (Heavier wire diameter)

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UCL	13 mm	3-0	VIOLET	70		36	V113H
		0	VIOLET	70		36	J1150H
 UCLX	16 mm	3-0	VIOLET	70		36	V13H
		2-0	VIOLET	70		36	V14H
		0	VIOLET	70		36	V15H
		1	VIOLET	70		36	V16H
 CT-3	22 mm	3-0	VIOLET	70		36	V327H
		2-0	VIOLET	70		36	V328H
		0	VIOLET	70		36	V329H
 CT-2	26 mm	3-0	VIOLET	70		36	V332H
		2-0	VIOLET	70		36	V333H
		2-0	UNDYED	70		36	V269H
		2-0	VIOLET	4x45 C/R		24	V7642E
		0	VIOLET	70		36	V334H
		0	UNDYED	70		36	V270H
		0	VIOLET	8x45 C/R		24	V744E
		0	VIOLET	4x70 C/R		24	4PK2575YE
		1	VIOLET	70		36	V335H
2	VIOLET	70		36	V336H		
 MO-6	27 mm	1	VIOLET	90	Heavy needle	36	V237H
 MO-5	31 mm	0	VIOLET	75	Heavy needle	12	W9215
		1	VIOLET	75	Heavy needle	12	W9216
 CT-1	36 mm	3-0	VIOLET	70		36	V338H
		3-0	UNDYED	90		36	V944H
		2-0	VIOLET	70		36	V339H
		2-0	VIOLET	90		36	V345H
		2-0	UNDYED	90		36	V945H
		2-0	UNDYED	8x70 C/R	Non-stocked item	12	JJ42D
		0	VIOLET	70		36	V340H
		0	VIOLET	90		36	V346H
		0	UNDYED	90		36	V946H
		0	VIOLET	8x45 C/R		12	J740D
		1	VIOLET	70		36	V341H
		1	VIOLET	90		36	V347H
		1	VIOLET	4x70 C/R		12	V736G
		2	VIOLET	90		36	V348H
2	UNDYED	90		36	V948H		
 MO-4	36 mm	0	VIOLET	70	Heavy needle	36	J436H
		1	VIOLET	90	Heavy needle	36	V255H





VICRYL* (Polyglactin 910) Braided or Monofilament Suture

1/2 Circle Taper Point Needles (Heavier wire diameter)


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CT	40 mm	3-0	VIOLET	70		36	V350H
		2-0	VIOLET	70		36	V351H
		2-0	VIOLET	8x45 C/R		12	8PM2544YG
		0	VIOLET	70		36	V352H
		0	VIOLET	90		36	V358H
		1	VIOLET	70		36	V353H
		1	VIOLET	90		36	V359H
		1	VIOLET	8x45 C/R		12	J753D
		2	VIOLET	90		36	V360H
	2	UNDYED	90		36	V960H	
 MO-45	45 mm	0	VIOLET	75	Heavy needle	12	W9244
		1	VIOLET	75	Heavy needle	12	W9245
		2	VIOLET	75	Heavy needle	12	W9246
 CTX	48 mm	2-0	VIOLET	70		36	V363H
		0	VIOLET	70		36	V364H
		0	VIOLET	90		36	V370H
		0	UNDYED	90		36	V170H
		0	VIOLET	8x45 C/R		12	J764D
		1	VIOLET	70		36	V365H
		1	VIOLET	90		36	V371H
		1	UNDYED	90		36	V366H
		1	VIOLET	8x45 C/R		12	J765D
		2	VIOLET	90		36	V372H
	2	UNDYED	90		36	V367H	
 TP-1	65 mm	2	VIOLET	90		12	V585G
 XLH	70 mm	2-0	VIOLET	70		12	V581G
		0	VIOLET	70		12	V582G

VICRYL* (Polyglactin 910) Braided or Monofilament Suture





1/2 Circle ETHIGUARD* Blunt Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 MHD-1	31 mm	3-0	VIOLET	75		12	W9983
		2-0	VIOLET	75		12	W9984
 CTD-1	36 mm	2-0	VIOLET	75		12	W9985
		2-0	VIOLET	75	Heavy Needle	12	W9993
		0	VIOLET	90	Heavy Needle	12	W9994
		1	VIOLET	90	Heavy Needle	12	W9995
 CTD	40 mm	0	VIOLET	75		12	W9988
		0	VIOLET	90	Heavy Needle	12	W9996
		1	VIOLET	90	Heavy Needle	12	W9997
 CTXD	48 mm	0	VIOLET	90		12	W9990
		1	VIOLET	90		12	W9991

3/8 Circle Blunt Point Needles



Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 BP-5	65 mm	1	VIOLET	100		12	W9391

5/8 Circle Taper Point Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UR-6	26 mm	3-0	VIOLET	70		36	J601H
		2-0	VIOLET	70		36	J602H
		0	VIOLET	70		36	V603H
		1	VIOLET	70		36	J2575H
 UR-6	26 mm	2-0	VIOLET	70		36	VCC2573H
 UR-5	36 mm	2-0	VIOLET	70		36	V375H
		0	VIOLET	70		36	V376H
		1	VIOLET	70		36	J377H
 UR-4	40 mm	0	VIOLET	70		36	J381H

VICRYL* (Polyglactin 910) Braided or Monofilament Suture






Fish Hook Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PB-30	30 mm	2-0	VIOLET	70	Taper point	36	DF2570H
 V-38	31 mm	1	VIOLET	75	TAPERCUT*	12	W9335

Straight Taper Point Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 ST-4	19 mm	3-0	VIOLET	70		36	BA2574H
 ST-4	19 mm	4-0	VIOLET	35		36	V620H
 ST	48 mm	2-0	VIOLET	90		36	V399H

1/2 Circle TAPERCUT* Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 V-5	17 mm	4-0	VIOLET	70		12	V994G
		4-0	VIOLET	70		36	V994H
		4-0	UNDYED	70		36	V230H
		3-0	VIOLET	70		36	V276H
 V-6	22 mm	4-0	UNDYED	70		36	V231H
 V-7	26 mm	3-0	VIOLET	70		12	V998G
		3-0	VIOLET	70		36	V998H
		2-0	VIOLET	70		36	V999H
		2-0	VIOLET	6x45 C/R		24	6OB2545YE
		0	VIOLET	70		36	V279H
 V-30	31 mm	2-0	VIOLET	75		12	W9360
		0	VIOLET	75		12	W9361
		1	VIOLET	75		12	W9362
 V-34	36 mm	3-0	VIOLET	90		36	V516H
		2-0	VIOLET	90		36	V517H
		0	VIOLET	90		36	V518H
		1	VIOLET	90		36	V519H

VICRYL* (Polyglactin 910) Braided or Monofilament Suture



1/2 Circle TAPERCUT* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 V-37	40 MM	2-0	UNDYED	75	Slimmed	12	W9900
		2-0	UNDYED	75		12	W9904
		2-0	VIOLET	90		12	W9466
		0	VIOLET	90		36	VS4000H
 V-39	45 mm	0	VIOLET	75		12	W9376
		1	VIOLET	75		12	W9377
 V-40	48 mm	2-0	VIOLET	90		36	VS5020H
		1	VIOLET	150	Loop	24	V829E
		6	VIOLET	75		12	V140G

3/8 Circle TAPERCUT* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CC-1	13 mm	6-0	UNDYED	45	MULTIPASS*	12	JV1727
 CC-1	13 mm	6-0	VIOLET	75	MULTIPASS*	36	JV1012
		5-0	VIOLET	75	MULTIPASS*	36	JV1013
 V-4	17 mm	5-0	UNDYED	75		12	W9824
		4-0	UNDYED	70		12	V240G
		3-0	UNDYED	75		12	W9826

5/8 Circle TAPERCUT* Needles







Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UV-1C	10 mm	6-0	VIOLET	30	ETHALLOY*	12	W9172
 UV-17	17 mm	4-0	UNDYED	75		12	W9820

VICRYL* (Polyglactin 910) Braided or Monofilament Suture

Compound Curve TAPERCUT* Needles



Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 DA	25 mm	4-0	VIOLET	45	VISI-BLACK*	12	18531G

1/2 Circle Reverse Cutting Needles









Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 X-1	22 mm	3-0	VIOLET	70		36	V460H
		3-0	UNDYED	70		36	ME2673H
		2-0	VIOLET	70		36	V461H
		2-0	UNDYED	75		36	JV459
 OS-4	22 mm	0	UNDYED	70		36	J694H
		2	VIOLET	75		12	W9287
 CP-2		2-0	VIOLET	75		36	JV969
		2-0	UNDYED	70		36	V869H
		0	VIOLET	70		36	MF2574H
		1	UNDYED	70		36	MF2675H
 CP-1		2-0	VIOLET	70		36	V466H
		2-0	VIOLET	90		36	V472H
		2-0	UNDYED	90		36	JV1472
		0	VIOLET	70		36	V467H
		0	VIOLET	90		36	V473H
		1	VIOLET	70		36	V468H
		1	VIOLET	90		36	V474H
		2	VIOLET	90		36	V482H
 CP	40 mm	2-0	VIOLET	90		36	JV484
		0	VIOLET	70		36	V479H
		0	VIOLET	90		36	V485H
		1	VIOLET	70		36	V480H
		1	VIOLET	90		36	V486H
 CPX	48 mm	0	VIOLET	90		36	V1057H
		1	VIOLET	90		36	V1058H
		2	VIOLET	90		36	V1059H

VICRYL* (Polyglactin 910) Braided or Monofilament Suture

1/2 Circle Reverse Cutting PRIME* Needles








Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 P-2	8 mm	5-0	UNDYED	45	ETHALLOY*	36	J503G
 PS-4	16 mm	5-0	UNDYED	45	Compound curve, ETHALLOY*	12	W9436
		4-0	UNDYED	45	MULTIPASS* ETHALLOY*	12	J507G
		4-0	UNDYED	45	Compound curve, ETHALLOY*	12	W9437

3/8 Circle Reverse Cutting PRIME* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 P-1	11 mm	6-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV489H
		6-0	UNDYED	45	ETHALLOY*	36	V489H
		5-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV490H
		5-0	UNDYED	45	ETHALLOY*	36	V490H
 P-3	13 mm	6-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	JV492H
		6-0	UNDYED	45		36	V492H
		5-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV493H
		5-0	UNDYED	45		36	V493H
		4-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV494H
		4-0	UNDYED	45		36	V494H
 PS-3	16 mm	5-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV500H
		4-0	UNDYED	45		24	W9506T
 PS-2	19 mm	5-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV495H
		4-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV496H
		4-0	UNDYED	45		24	W9515T
		3-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV497H
		3-0	UNDYED	45		36	V497H
 PS-2 S	19 mm	4-0	UNDYED	45	Slimmed, ETHALLOY*	36	V496H
 PS-1	24 mm	4-0	UNDYED	70		36	JE2672H
		3-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV683H
		3-0	UNDYED	45		36	V683H
		2-0	UNDYED	45	MULTIPASS* ETHALLOY*	36	MPV684H
		2-0	UNDYED	45		36	V684H
 PS	26 mm	3-0	UNDYED	70	MULTIPASS* ETHALLOY*	36	MPV9526H
		2-0	UNDYED	75		24	W9527T
 PSL	30 mm	2-0	VIOLET	70		36	V686H





VICRYL* (Polyglactin 910) Braided or Monofilament Suture

3/8 Circle Reverse Cutting Needles


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
	13 mm	6-0	VIOLET	45		36	J384H
		5-0	VIOLET	45		36	V385H
		4-0	VIOLET	45		36	V386H
	13 mm	4-0	VIOLET	45		12	JV559
		5-0	VIOLET	45		36	V388H
		4-0	VIOLET	75	CONV. CUTTING	12	W9443
	19 mm	5-0	UNDYED	45		36	V291H
		5-0	VIOLET	45		36	V391H
		5-0	UNDYED	70		36	V421H
		4-0	UNDYED	45		36	V292ZH
		4-0	VIOLET	45		12	V392ZG
		3-0	UNDYED	45		36	V293H
		3-0	VIOLET	45		12	V393G
		3-0	VIOLET	45		36	V393H
		3-0	UNDYED	70		36	V423H
		3-0	VIOLET	70		36	V398H
	19 mm	4-0	UNDYED	45	Slimmed	36	V292H
		4-0	VIOLET	45	Slimmed	36	V392H
		4-0	UNDYED	70	Slimmed	36	V422H
		4-0	VIOLET	70	Slimmed	36	V397H
	24 mm	4-0	UNDYED	70		36	V441H
		4-0	VIOLET	70		36	V451H
		3-0	UNDYED	70		36	V442H
		3-0	VIOLET	70		36	V452H
		2-0	UNDYED	70		36	V443H
		2-0	VIOLET	70		36	V453H
	30 mm	3-0	UNDYED	70		36	V588H
		3-0	VIOLET	70		36	V585H
		2-0	UNDYED	75		36	JV589
		2-0	VIOLET	70		36	V586H
		0	UNDYED	70		36	V520H
		0	VIOLET	70		36	V587H

VICRYL* (Polyglactin 910) Braided or Monofilament Suture



3/8 Circle Conventional Cutting PRIME* Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 PC-3	16 mm	5-0	UNDYED	45	MULTIPASS* ETHALLOY*	24	W9505T
		4-0	UNDYED	45	MULTIPASS* ETHALLOY*, non-stocked	36	MPV9506H
 PC-5	19 mm	4-0	UNDYED	75	ETHALLOY*	24	W9570T
		3-0	UNDYED	75	MULTIPASS* ETHALLOY*	24	W9571T
 PC-25	26 mm	4-0	UNDYED	70	MULTIPASS* ETHALLOY*, non-stocked	36	MPV9580H
		4-0	UNDYED	75		24	W9580T
		3-0	UNDYED	75		24	W9581T
 PCLX	36 mm	3-0	UNDYED	75	ETHALLOY*	24	W9531T

Straight Cutting Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 SC-1	13 mm	4-0	VIOLET	45	ETHALLOY*	36	V383H
 KS	60 mm	4-0	UNDYED	70		36	J662H
		3-0	VIOLET	70		36	V523H
		3-0	UNDYED	70		36	V663H
		2-0	VIOLET	75		12	W9418
		2-0	UNDYED	75		12	W9718

1/2 Circle MICROPOINT* Spatula

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 GS-14 GS-14	5,5 mm	10-0	VIOLET	30		12	V450G
		9-0	VIOLET	30		12	V449G
		8-0	VIOLET	20		12	V574G
		7-0	VIOLET	30		36	ENN2534H
 TG175-8	7 mm	7-0	VIOLET	45	ETHALLOY*	12	W9565

VICRYL* (Polyglactin 910) Braided or Monofilament Suture





3/8 Circle MICROPOINT* Spatula

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 GS-9	6,5 mm	8-0	VIOLET	30		12	V120G
		7-0	VIOLET	30		12	V119G
 GS-9 GS-9	6,5 mm	10-0	VIOLET	30	Monofilament	12	V440G
		9-0	VIOLET	30		12	V549G
		9-0	VIOLET	30	Monofilament	12	V439G
		8-0	VIOLET	20		12	V547G
		8-0	VIOLET	30		12	V548G
		7-0	VIOLET	45		12	V546G
 TG140-E TG140-E	6,6 mm	7-0	VIOLET	30		12	W9561

1/4 Circle MICROPOINT* Spatula Needles


Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 TG100-8	6,5 mm	6-0	VIOLET	45	ETHALLOY*	12	V544G

3/8 Circle MICROPOINT* Reverse Cutting Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 TG140-4	6,2 mm	10-0	VIOLET	30	Monofilament	12	EMM2830G
 G-6	8 mm	7-0	VIOLET	30		36	EA2531H
 G-6 G-6	8 mm	6-0	VIOLET	45	ETHALLOY*	12	V551G
 G-1 G-1	11 mm	6-0	VIOLET	45	ETHALLOY*	12	W9537
 G-3 G-3	13 mm	5-0	VIOLET	45	ETHALLOY*	12	V557G

VICRYL* (Polyglactin 910) Braided or Monofilament Suture


3/8 Circle CS ULTIMA* Spatula

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
CS-140-6 	6,2 mm 10-0	VIOLET	15		12	V429G


1/8 Circle CS ULTIMA* Spatula

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
CS44-6 	6,5 mm 10-0	VIOLET	10	Monofilament	12	V965




3/8 Circle Spatula Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
SPAT-10 	11 mm	5-0	UNDYED	45	ETHALLOY*	12	W9761
		4-0	UNDYED	45	ETHALLOY*	12	W9762

1/2 Circle Spatula Needles






Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
S-22 S-22 	8 mm	6-0	VIOLET	45		12	W9566
		6-0	UNDYED	45		12	W9756
		5-0	VIOLET	45		12	W9567

1/4 Circle Spatula Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
S-4 S-4 	8 mm 450 microns 4-0	VIOLET	45	ETHALLOY*	12	V572G	
S-29 S-29 	8 mm 225 microns 6-0	VIOLET	45		12	J555G	
S-14 S-14 	8 mm, 350 microns	6-0	VIOLET	45	ETHALLOY*	12	V570G
		6-0	VIOLET	45	VISI-BLACK* ETHALLOY*	12	V575G
		6-0	UNDYED	45	ETHALLOY*	12	V670G
		6-0	VIOLET	20	ETHALLOY*	12	W9562
		5-0	VIOLET	45		12	V571G
		5-0	UNDYED	45		12	V671G

VICRYL* (Polyglactin 910) Braided or Monofilament Suture

3/8 Circle Taper Point Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 BV-6	4,8 mm	10-0	VIOLET	10	Monofilament	12	V970G
 BV-4	4,9 mm	10-0	VIOLET	15	Monofilament	12	V971G
 BV130-5	6,5 mm	9-0	VIOLET	15	Monofilament, MULTIPASS*	12	V239G
		8-0	VIOLET	20	MULTIPASS*	12	V400G
 BV-1	9,3 mm	8-0	VIOLET	30		12	W9577
		6-0	VIOLET	30		12	W9575
 BB	17 mm	4-0	VIOLET	75	VISI-BLACK*	36	JV2024
		3-0	VIOLET	75	VISI-BLACK*	36	JV2025

VICRYL* (Polyglactin 910) Braided or Monofilament Suture

Ligatures

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
LIGAPAK*	4-0	VIOLET	250	Reel	12	V1204G
	3-0	VIOLET	250	Reel	12	V1205G
	2-0	VIOLET	250	Reel	12	V1206G
	0	VIOLET	250	Reel	12	V1207G
	1	VIOLET	250	Reel	12	V1208G
	3-0	VIOLET	125	Reel	12	V205G
	2-0	VIOLET	125	Reel	12	V206G
SUTUPAK*	3-0	VIOLET	140		36	V614H
	2-0	VIOLET	140		36	V615H
	0	VIOLET	140		36	V616H
	1	VIOLET	140		36	V617H
	2	VIOLET	140		36	V618H
SUTUPAK*	4-0	VIOLET	2x70		24	V623E
	3-0	VIOLET	2x70		24	V624E
	2-0	VIOLET	2x70		24	V625E
	0	VIOLET	2x70		24	V626E
	1	VIOLET	2x70		36	V627H
	2	VIOLET	2x70		36	V628H
SUTUPAK*	5-0	VIOLET	3x45		36	V632H
	4-0	VIOLET	3x45		36	V633H
	4-0	UNDYED	3x45		36	V643H
	3-0	VIOLET	3x45		36	V634H
	3-0	UNDYED	3x45		36	V644H
	2-0	VIOLET	3x45		36	V635H
	2-0	UNDYED	3x45		36	V645H
	0	VIOLET	3x45		24	V636E
	1	VIOLET	3x45		24	V637E
	2	VIOLET	3x45		36	V638H

VICRYL* (Polyglactin 910) Braided or Monofilament Suture

Ligatures

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
SUTUPAK*	4-0	VIOLET	5x70		24	V1214E
	3-0	VIOLET	5x70		24	V1215E
	3-0	UNDYED	5x70		24	V963E
	2-0	VIOLET	5x70		24	V1216E
	2-0	UNDYED	5x70		25	V964E
	0	VIOLET	5x70		24	V1217E
	0	UNDYED	5x70		24	V965E
	1	VIOLET	5x70		36	V1218H
	2	VIOLET	5x70		36	V1219H
SUTUPAK*	5-0	VIOLET	6x45		36	V1223H
	4-0	VIOLET	6x45		36	V1224H
	3-0	VIOLET	6x45		36	V1225H
	3-0	UNDYED	6x45		36	V978H
	2-0	VIOLET	6x45		36	V1226H
	0	VIOLET	6x45		24	V1227E
	1	VIOLET	6x45		24	V1228E
	2	VIOLET	6x45		36	V1229H
SUTUPAK*	4-0	VIOLET	12x45		36	V903H
	3-0	VIOLET	12x45		36	V904H
	3-0	UNDYED	12x45		36	V910H
	2-0	VIOLET	12x45		24	V905E
	2-0	UNDYED	12x45		24	V911E
	0	VIOLET	12x45		24	V906E
	0	UNDYED	12x45		24	V912E
	1	VIOLET	12x45		24	V907E