

- NEW PRODUCT INFORMATION SHEET -

Specifications subject to change
© OMRON HEALTHCARE EUROPE B.V.

Name: Flex Temp Smart

Thermometer Model: MC-343F-E



Product Features:

- Digital Fever Thermometer for Oral, Axillary and Rectal use
- Fast 10 seconds measurement (Rectal: approx.10 seconds, Oral: approx. 20 seconds, Armpit: approx. 25 seconds)
- Flexible tip for highest comfort during measurement, especially for children/babies
- Protective case to highlight the value of the device
- Beeper
- Celsius/Fahrenheit indication
- 3-digit display (1 decimal place)
- Water-resistant
- Memory (last measurement value)
- Exchangeable battery
- Unique and attractive design
- 15 language Instruction Manual

Product Position: Premium product for design oriented segment of the market

The **Flex Temp Smart** is a thermometer following an award winning style of OMRON's other thermometer styles. The special flexible tip offers a soft comfort that will enable the patient to use the thermometer safely and conveniently during all measurement methods. Being water-resistant and having an exchangeable battery are natural features for this product. The easy-to-use placement of the On/Off-button and the excellent design reinforces the recognition as a high-quality OMRON-thermometer.

The packaging will be a "blister" package where you can clearly see the product. Very clean, stylish packaging to compliment the OMRON thermometer family design.

Package example:

Front Design:



Back Design:



Complete Package:



Package will have "Rectal: approx. 10 seconds" on the back as the explanation of "10 Seconds Technology".

Instruction Manual Layout example:



Pictures available on Sample Kit CD:



Languages for both the Instruction Manual and Packaging:

15 languages

English, German, French, Italian, Spanish, Portuguese, Danish, Finnish, Norwegian, Swedish, Dutch, Greek, Turkish, Russian, Arabic

OMRON Flex Temp Smart product specifications:

Product Type and Number	Digital Thermometer MC-343F-E
Product Name	Flex Temp Smart
Power Supply	1.5V DC, 1LR41 Alkaline-Manganese Button Battery
Power Consumption	0.1mW
Battery Life	Approx. 2 years or more (3 times per day)
Sensing Unit	Thermistor
Measurement Method	Actual measurement (non-predictive)
Measurement Range	32.0°C to 42.0°C (89.6°F to 107.6°F) Lo < 32.0°C (89.6°F) Hi > 42.0°C (107.6°F)
Temperature Scales	C and F (switchable)
Display Resolution	0.1°C or 0.1°F
Accuracy	±0.1°C between 32.0°C to 42.0°C ±0.2°F between 89.6°F to 107.6°F
Measurement Position	Oral, Axillary (armpit), and Rectal
Auto Shut-Off Time	Thirty minutes after measurement reading or three minutes when turned on but not used
Memory Recall	Last temperature
Operating Temperature / Humidity	+10.0°C to +40.0°C (+50.0°F to +104.0°F) 30% to 85% RH
Storage Temperature / Humidity	-20.0°C to +60.0°C (-4.0°F to + 140.0°F) 10% to 95% RH
Dimensions	19.4mm (W) x 132.5mm (L) x 10.0mm (D)
Weight	Approx. 11 grams (including battery)
Package Content	Test Battery (LR41), Storage Case, Instruction Sheet
Warranty	Three Years
EAN code	401567210491 4

OMRON Flex Temp Smart packaging specifications:

Carton	QTY	Approx. Weight	Approx. Dimensions L*W*D (mm)
Unit (without pkg)	1	17.2 g	132.5 x 19.4 x 10
Individual Package	1	70.7 g	180.9 x 76.7 x 23.9
Inner carton	10	840.6 g	277 x 87 x 183
Master carton	100	10.6 kg	570 x 455 x 214

Takefumi Nakanishi
Product Technical Manager
OMRON HEALTHCARE EUROPE B.V.