



OMRON M6 Comfort comes with the Intelli Wrap Cuff technology, which supports 360 degrees accuracy, enabling accurate results in any position around the upper arm¹.

- Stroke prevention – Helps to detect stroke risk¹
- Clinically validated
- Intellisense Technology – measurement without unpleasant feeling of pressure
- Intelli Wrap Cuff Technology – 360 accuracy. Accurate results in any position around the upper arm.²
- Blood Pressure Level Indicator – Indicates if result is above recommended level
- Cuff Wrap Guide – Supports accurate measurement
- Irregular Heartbeat Detection
- Averaging Function – Calculates the average of most recent values
- Morning Averaging Function – Calculates the average of up to three measurements taken within 10 minutes of the first measurement of the morning (4:00- 11:59)
- Evening Averaging Function – Calculates the average of up to three measurements taken within 10 minutes of the first measurement of the evening (19:00- 1:59)
- Advanced Averaging Function – Weekly morning and evening averages
- Body Movement Detection
- Morning Hypertension Tracker – Indicates if morning average is above recommended level
- Memory Capacity – 2 x 100
- Storage bag included