

SILON monofilament

SILON monofilament is a non-absorbable, synthetic suture made of polyamide, which has a smooth non-porous surface.

The thread is extremely useful for skin sutures and has a high tensile strength in knots, and knots tie down easily. The thread is available in blue and black.

Silon monofilament can cause a slight initial inflammatory tissue reaction.



Indication: Skin Closure (e.g. intracutaneous, subcutaneous), Plastic, General, Reconstructive and Ophthalmic surgery.

Contraindication: Silon monofilament should not be used when long-term mechanical support is required.

Chemical Composition: Polyamide 6

Sizes: 10/0 to 2 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Blue and Black

SILON monofilament

UK03
160



Non-absorbable - Silon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 DS 15	0.7	6/0	1 x 45 cm	Blue	SM 2260	
	1	5/0	1 x 45 cm	Blue	SM 2261	
	1.5	4/0	1 x 45 cm	Blue	SM 2262	
	2	3/0	1 x 45 cm	Blue	SM 2263	
	1.5	4/0	1 x 75 cm	Blue	SM 2362	
 DS 16	0.7	6/0	1 x 45 cm	Blue	SM 2460	
	0.7	6/0	1 x 60 cm	Blue	SM 2459	
 DS 19	0.7	6/0	1 x 45 cm	Blue	4332	
	1	5/0	1 x 45 cm	Blue	SM 2266	
	1.5	4/0	1 x 45 cm	Blue	SM 2367	
	2	3/0	1 x 45 cm	Blue	SM 2368	
	3	2/0	1 x 45 cm	Blue	SM 2369	
	1	5/0	1 x 75 cm	Blue	SM 2366	
	1.5	4/0	1 x 75 cm	Blue	SM 2267	
	2	3/0	1 x 75 cm	Blue	SM 2268	
	3	2/0	1 x 75 cm	Blue	SM 2269	
	1	5/0	1 x 90 cm	Blue	SM 2466	
 DS 20	1.5	4/0	1 x 75 cm	Blue	SM 2471	
 DS 24	1.5	4/0	1 x 75 cm	Blue	SM 3277	
	2	3/0	1 x 75 cm	Blue	SM 3276	
	3	2/0	1 x 75 cm	Blue	SM 3275	
 DS 25	1.5	4/0	1 x 45 cm	Blue	SM 2372	
	2	3/0	1 x 45 cm	Blue	SM 2373	
	3	2/0	1 x 45 cm	Blue	SM 2374	
	1	5/0	1 x 75 cm	Blue	SM 2375	
	1.5	4/0	1 x 75 cm	Blue	SM 2273	
	2	3/0	1 x 75 cm	Blue	SM 2270	
	3	2/0	1 x 75 cm	Blue	SM 2271	
	3.5	0	1 x 75 cm	Blue	SM 2272	
	4	1	1 x 75 cm	Blue	SM 2274	
	2	3/0	1 x 90 cm	Blue	SM 3272	
	3	2/0	1 x 90 cm	Blue	SM 3271	