

Gauze swabs, 5-pack, sterile

REF: 222294, 222295, 222296, 2222765

SELEFA® Gauze swabs, sterile are used for non-invasive medical procedures, such as wound paddings, absorbing blood and exudates, skin and wound cleansing or as secondary dressing.

Key Features:

- ▶ 100% cotton gauze (17 threads per cm²)
- ▶ Peroxide bleached (no chlorine)
- ▶ 12-ply folded
- ▶ Folded in cut edges to prevent loose threads
- ▶ Sterilized by ethylene oxide gas (EO)
- ▶ 5 swabs per pouch



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
222294	Gauze swabs, 5-pack, 12-ply, sterile	10 x 10 cm	80/1280
222295	Gauze swabs, 5-pack, 12-ply, sterile	7,5 x 7,5 cm	80/1600
222296	Gauze swabs, 5-pack, 12-ply, sterile	5 x 5 cm	80/1600
2222765	Gauze swabs, 5-pack, 12-ply, sterile	10 x 20 cm	40/720

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices
EO sterilized according to valid version of ISO 11135.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5 and EN ISO 10993-10

Compliance with legal requirements: Medical Devices Directive 93/42/EEC, Class Is

Compliance with product standards: EN 14079

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1 and EN 1041

Material

Product Part	Raw Material & Characteristic
Raw material	100 % cotton gauze, 17 threads/cm ² , white

Warnings and Recommendations

Warnings:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical paper / Plastic film PA/PE	222294: 75 x 90
		222295: 100 x 114
		222296: 133 x 135
		2222765: 133x241

Dispenser Pack	Carton, ECF or TCF quality	222294: 90 x 78 x 145 222295: 115 x 95 x 145 222296: 135 x 130 x 145 2222765: 250 x 130 x 90
Transport Carton	Carton, ECF or TCF quality	222294: 400 x 190 x 300 222295: 485 x 240 x 300 222296: 530 x 280 x 300 2222765: 400 x 280 x 510
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
222294	 6 438129 151442	 6 438129 351446	 6 438129 551440
222295	 6 438129 151459	 6 438129 351453	 6 438129 551457
222296	 6 438129 151466	 6 438129 351460	 6 438129 551464
2222765	 6 438129 151046	 6 438129 351040	 6 438129 551044

ONEMED

Manufactured in China for
OneMed Group Oy
 Address: Metsäläntie 20
 Postal Address: P.O. Box 10
 FI-00321 Helsinki, Finland
 Phone +358 20 786 6800
 www.onemed.com
 Domicile Helsinki
 Business ID 2039640-1