

LIMITING EXPOSURE AT EVERY STEP

Nonin WristOx₂[®] 3150

Wrist-Worn Pulse Oximeter



6000CA-WO2 Disposable
Cloth Pulse Oximetry Sensor



Tamper-proof wristband



WristOx₂[®] 3150
Pulse Oximeter*



WristOx₂[®] 3150 System

Introducing Single-Patient Accessories for the WristOx₂[®] 3150

Preventing infection during a pandemic means you need safety at every step. Nonin single-patient sensors and wrist straps reduce touchpoints to help limit exposure while capturing reliable, accurate patient data. **Safety first. Reliability always.**

*WristOx₂ 3150 is sold in two configurations based on desired data download/accessibility. The WristOx₂ 3150 USB support data download through USB connection and the WristOx₂ 3150 with BLE supports data transmitted via Bluetooth Low Energy.



Single Patient Use, Multiple Infection Prevention Benefits



Reduces spread of infection from patient to patient and patient to clinician



PureLight LEDs eliminate variation in readings for consistency across studies and patients



Maximizes workflow and time-savings by reducing number of components requiring cleaning between studies



Works on all skin pigmentations for quality monitoring across all patient populations



Eliminates need for cleaning accessories and potential variability in cleaning practices and efficacy



Every Nonin sensor is functionally tested for guaranteed quality right out of the box

Designed With Infection Prevention in Mind

Eliminates need for cleaning accessories between use

PureLight LEDs eliminate variation in readings

Short cable length designed for wrist-to-finger fit

Breathable, flexible sensor cloth material

Tamper-proof wristband closure

WristOx²® 3150
Pulse Oximeter



6000CA-WO2 Disposable Cloth
Pulse Oximetry Sensor



Tamper-proof
wristband

Nonin.com

United States
info@nonin.com
+1-763-553-9968

International
infointl@nonin.com
+31 (0)13-45 87 130

