

Patient and standing up hoist



www.rebotec.de

Arnold 150

Patient hoist



ref. no. 320.00.10

Equipment:

- electric lifting arm
- 24 Volt motor / dry battery
- cable remote control
- emergency stop
- emergency lowering
- lowering to the ground is possible
- low construction of the base frame
= good underride possibility
- undercarriage is expandable
by means of pedal
- swivelling castors, rear castors with brake
- assembled in a few minutes
- approved up to 150 kg body weight

ref. no. 320.00.10

Arnold 150 with scale PW600

as a set

ref.no. 320.06.10


Accessories:

Four-point-sling-hanger, wide
for Arnold 150

- 56 cm x 33 cm x 14,5 cm
ref. no. 420.10.42



ref.-no 420.10.42



Type:	Arnold 150
max. body weight:	150 kg
lift height (min):	55 cm
lift height (max):	174 cm
underpin height:	13 cm
weight:	43 kg
ref. no.:	320.00.10



For patient hoists


1. Universal sling:

- high back with removable PE-strips
- upholstered leg parts
- water resistant, washable
- approved up to 250 kg body weight

size S ref. no. 321.10.10

size M/L ref. no. 321.10.20

size XL/XXL ref. no. 321.10.30



2. Universal sling with integrated head rest

- approved up to 300 kg body weight

size M/L ref. no. 321.40.20

www.rebotec.de

Distributed by:

REBOTEC® Rehabilitationsmittel GmbH

D-49610 Quakenbrück · Artlandstraße 57 – 59 · Fon +49 (0) 5431 9416-0 · Fax +49 (0) 5431 9416-66

www.rebotec.de · info@rebotec.de