

SPECIFICATION

Pocket Breezer



The foldable pocket mask ensures the safety of the lifesaver giving artificial respiration.

The mouthpiece has a bacterial filter and a one-way valve for the safety of the patient. The valve is positioned between the mouthpiece and air-cushioned mask to reduce the risk of infection.

The close-fitting mask possesses an attachment strap and folds into itself for easy storage. Additionally, it has got a bacterial filter for the patients safety.

Comes in a practical Lifeguard pocket-box.

Comes in a practical Lifeguard pocket-box.

SPECIFICATION

Pocket Breezer

Type	PZN	REF	Unit
Pocket Breezer, complete		N4 2050	1 pc.

