

Non-woven swabs

Compressas em não tecido

Compresses en non tissé

Vliesstoff-Kompressen

Nonwoven kompres

Compressa in T.N.T.

Compresa de tela sin tejer

Non-woven kompress



Mais imagens



Sterile non-woven swabs

- For wound dressing and cleansing. For treatment of lightly or moderately exuding wounds.
- 70% viscose, 30% polyester.
- Approx. weight: 30, 40 and 50g/sq.m.
- Some models with X-Ray detectable thread with approx. 60% barium sulfate contents.
- Some models available with Y shape cut.
- Peel-open pack + shelf box + transport carton.



Code	Size	Plies	g/m2	X-Ray	Y shape cut	Pcs./Pack	Pcs./Sh.box	Pcs./Tr. carton
482-001	5 x 5cm	4	30			1	160	4320
482-002	5 x 5cm	4	30			2	240	6480
482-003	5 x 5cm	4	30			3	240	6480
482-004	5 x 5cm	4	30			5	250	6750
482-005	5 x 5cm	4	30			10	280	7560
482-006	7,5 x 7,5cm	4	30			1	160	2880
482-007	7,5 x 7,5cm	4	30			2	240	4320
482-008	7,5 x 7,5cm	4	30			3	240	4320
482-009	7,5 x 7,5cm	4	30			5	250	4500
482-010	7,5 x 7,5cm	4	30			10	280	5040
482-016	10 x 10cm	4	30			1	160	1600
482-017	10 x 10cm	4	30			2	240	2400
482-018	10 x 10cm	4	30			3	240	2400
482-019	10 x 10cm	4	30			5	250	2500
482-020	10 x 10cm	4	30			10	280	2800
482-021	10 x 20cm	4	30			1	160	960
482-022	10 x 20cm	4	30			2	240	1440
482-023	10 x 20cm	4	30			3	240	1440
482-024	10 x 20cm	4	30			5	250	1500
482-025	10 x 20cm	4	30			10	280	1680
482-026	5 x 5cm	4	40			1	150	4050
482-027	5 x 5cm	4	40			2	180	4860
482-029	5 x 5cm	4	40			5	225	6075
482-031	7,5 x 7,5cm	4	40			1	150	2700
482-032	7,5 x 7,5cm	4	40			2	180	3240
482-034	7,5 x 7,5cm	4	40			5	225	4050
482-041	10 x 10cm	4	40			1	150	1500
482-042	10 x 10cm	4	40			2	180	1800
482-044	10 x 10cm	4	40			5	225	2250
482-045	10 x 10cm	4	40			10	220	2200
482-046	10 x 20cm	4	40			1	150	900
482-047	10 x 20cm	4	40			2	180	1080
482-049	10 x 20cm	4	40			5	225	1350
482-144	10 x 10cm	4	50	X		5	200	2000
482-149	10 x 20cm	4	50	X		5	200	1200
482-162	15 x 20cm	4	30	X		5	400	2400
482-163	15 x 20cm	4	30			5	400	2400
482-211	7,5 x 7,5cm	4	30		X	2	160	2880
482-213	10 x 10cm	4	30		X	2	160	1600
482-226	5 x 5cm	4	40		X	5	200	5400
482-229	10 x 10cm	4	40		X	5	200	2000