




evercare® Adhesive aperture drape, transparent, 80 x 130 cm

REF: 1516

evercare® Adhesive aperture drape, sterile with adhesive aperture provides easy and quick draping for a variety of procedures.

Key Features

- Easy and quick draping
- Transparent
- Complies to EN 13795-1 High performance critical product area requirements



80 x 130 cm
Ø 15 x 6 cm


REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
1516	Adhesive aperture drape, sterile, Transparent	80 x 130 cm, Ø 15 x 6 cm	35 / 140

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-01 EN ISO 10993-05 EN ISO 10993-07 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN13795-1
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Cover material	PE, transparent

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, Transparent	150 x 290
Dispenser Pack	TCF or ECF quality	265 x 143 x 234
Transport Carton	TCF or ECF quality	545 x 310 x 246
Packing Tape	Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
1516	6438129118605	6438129318609	6438129518603

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.