

evercare® Elastic support bandage Ideal

REF: 4240-01, 4241-01, 4242-01, 4243-01, 4244-01, 4245-01, 4246, 4247, 4248, 4249, 4250

evercare® Elastic support bandage Ideal is stretchable with soft and porous structure, have high permeability and is comfortable to wear. Thanks to woven-edge weaving technique and careful selection of raw material, the bandages may be used for support, fixation of most dressings, especially for joints and round body parts, and small medical devices to the patient's body. Package designed for easy product identification.

Key Features

- ▶ Fixation of wound dressings, especially on joints and round body parts
- ▶ For support, e.g. after plaster removal and odema treatment
- ▶ Resistant to ointments
- ▶ Bandage with woven edges
- ▶ Soft and comfortable
- ▶ 95% stretchability
- ▶ 100 % cotton




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
4240-01	Elastic support bandage Ideal, cotton, non-sterile	4 cm x 4 m	10 / 480
4241-01	Elastic support bandage Ideal, cotton, non-sterile	6 cm x 4 m	10 / 400
4242-01	Elastic support bandage Ideal, cotton, non-sterile	8 cm x 4 m	10 / 320
4243-01	Elastic support bandage Ideal, cotton, non-sterile	10 cm x 4 m	10 / 240
4244-01	Elastic support bandage Ideal, cotton, non-sterile	12 cm x 4 m	10 / 240
4245-01	Elastic support bandage Ideal, cotton, non-sterile	15 cm x 4 m	10 / 160
4246	Elastic support bandage Ideal, cotton, non-sterile	6 cm x 4 m	10 / 320
4247	Elastic support bandage Ideal, cotton, non-sterile	8 cm x 4 m	10 / 240
4248	Elastic support bandage Ideal, cotton, non-sterile	10 cm x 4 m	10 / 200
4249	Elastic support bandage Ideal, cotton, non-sterile	12 cm x 4 m	10 / 160
4250	Elastic support bandage Ideal, cotton, non-sterile	15 cm x 4 m	10 / 120

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	DIN 61632, DIN 61634
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Bandage	100% cotton, natural white, Unstretched length 2 m, Stretched lenght 4 m, Elongation approx. 90%, Elasticity approx. 80%

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	4240-01: 235 x 95 x 45 4241-01: 235 x 95 x 65 4242-01: 235 x 95 x 85 4243-01: 235 x 95 x 105 4244-01: 235 x 95 x 125 4245-01: 235 x 95 x 155 4246: 240 x 100 x 65 4247: 240 x 100 x 85 4248: 240 x 100 x 85 4249: 240 x 100 x 125 4250: 240 x 100 x 125
Transport Carton	Carton, TCF or ECF quality	4240-01: 490 x 310 x 410 4241-01: 490 x 360 x 410 4242-01: 490 x 370 x 410 4243-01: 490 x 340 x 410 4244-01: 490 x 395 x 410 4245-01: 490 x 340 x 410 4246: 490 x 210 x 540 4247: 490 x 210 x 530 4248: 490 x 210 x 530 4249: 490 x 210 x 520 4250: 490 x 210 x 520
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4240-01	N/A	06438129321784	06438129521788
4241-01	N/A	06438129321791	06438129521795
4242-01	N/A	06438129321807	06438129521801
4243-01	N/A	06438129321814	06438129521818
4244-01	N/A	06438129321821	06438129521825
4245-01	N/A	06438129321838	06438129521832
4246	06438129114980	06438129314984	06438129514988
4247	06438129114997	06438129314991	06438129514995
4248	06438129115000	06438129315004	06438129515008
4249	06438129115017	06438129315011	06438129515015
4250	06438129115024	06438129315028	06438129515022

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.