

Product overview Qufora urine bags

Qufora Bed bags

REF	Product spec	Sterile	Retail	No. Per shipper
16823201	Turp bag 2000 ml, 125 cm inlet, sample port, drip chamber, big lever, int hanger	S	–	10
16923201	UriMeter 2600 ml, 130 inlet tube, sample port, drip chamber, big lever, int hanger	S	–	10
16923301	Turp bag 5000 ml, 125 cm inlet, sample port, drib chamber, big lever, int hanger	S	–	10
14823202	Bed Bag 2000 ml (U4), 120 cm inlet tube, big lever, air filter, sample port	S	–	25
14853202	Bed Bag 2000 ml (Q4), 120 cm inlet tube, free-flow, sample port	S	–	25
14851201	Bed Bag 2000 ml (S4), 130 cm inlet tube, big lever, air filter, sample port	S	–	25
13853101	Bed Bag 2000 ml, 90 cm inlet tube, free-flow	N	25	100
13853401	Bed Bag 2000 ml, 110 cm inlet tube, free-flow	N	25	100
13853402	Bed Bag 2000 ml, 110 cm wide inlet tube, free-flow	N	–	100
13850301	Bed Bag 2000 ml, 130 cm inlet tube, free-flow	N	25	100
13753101	Bed Bag 1500 ml, 90 cm inlet tube, free-flow	N	25	100
13303802	Bed Bag 750 ml, 10 cm inlet tube, no outlet	N	25	100
13803112	Bed Bag 2000 ml, 90 cm inlet tube, no outlet	N	25	100
13803412	Bed Bag 2000 ml, 110 cm inlet tube, no outlet	N	25	100
13703802	Bed Bag 1500 ml, 10 cm inlet tube, no outlet	N	25	100
13703902	Bed Bag 1500 ml, 35 cm inlet tube, no outlet	N	25	100
13703112	Bed Bag 1500 ml, 90 cm inlet tube, no outlet	N	25	100
13703412	Bed Bag 1500 ml, 110 cm inlet tube, no outlet	N	25	100

Qufora Leg bags - MULTI

REF	Product spec	Sterile	Retail	No. Per shipper
11061901	Leg Bag 350 ml, 30 cm inlet tube, extra connector, Q-flow valve, soft backing	N	10	100
11121501	Leg Bag 500 ml, multi, 50 cm adj inlet tube, T-valve, soft backing	N	10	100
11151901	Leg Bag 500 ml, multi, 30 cm inlet tube, extra connector, T-valve, soft backing	N	10	100
11161902	Leg Bag 500 ml, multi 30 cm inlet tube, extra connector, Q-flow, soft backing	N	10	100
11161502	Leg Bag 500 ml, multi, 50 cm adj inlet tube, Q-flow valve, soft backing	N	10	100
11321501	Leg Bag 750 ml, multi, 50 cm adj inlet tube, T-valve, soft backing	N	10	100
11361501	Leg Bag 750 ml, multi, 50 cm adj inlet tube, Q-flow valve, soft backing	N	10	100
12061901	Leg Bag 350 ml, 30 cm inlet tube, extra connector, Q-flow valve, soft backing	S	10	100
12121501	Leg Bag 500 ml, multi, 50 cm inlet tube, fixed connector T-valve	S	10	100
12121601	Leg Bag 500 ml, multi, 25 cm inlet tube, fixed connector T-valve	S	10	100
12121701	Leg Bag 500 ml, multi, 2 cm inlet tube, fixed connector T-valve	S	10	100
12161502	Leg Bag 500 ml, multi, 50 cm inlet tube, fixed connector + spare connector Q-flow valve	S	10	100
12161602	Leg Bag 500 ml, multi, 25 cm inlet tube, fixed connector + spare connector Q-flow valve	S	10	100
12161802	Leg Bag 500 ml, multi, 10 cm straight inlet tube, fixed connector, sample port, Q-flow valve	S	10	100
12161901	Leg Bag 500 ml, multi, 30 cm inlet tube, extra connector, Q-flow valve, soft backing	S	10	100
12321501	Leg Bag 750 ml, multi, 50 cm tube, fixed connector, T-valve	S	10	100
12321601	Leg Bag 750 ml, multi, 25 cm tube, fixed connector, T-valve	S	10	100
12361502	Leg Bag 750 ml, multi, 50 cm tube, fixed connector + spare connector, Q-flow valve	S	10	100

Qufora Leg bags - STANDARD

REF	Product spec	Sterile	Retail	No. Per shipper
11223802	Leg Bag 600 ml, standard, 10 cm angled inlet, small lever	N	25	100
11221902	Leg Bag 600 ml, standard, 30 cm inlet, small lever, soft backing	N	25	100
11221502	Leg Bag 600ml, standard, 50 cm adj tube, small lever, soft backing	N	25	100
11259802	Leg Bag 600 ml, standard, 10 cm angled inlet, T-valve	N	25	100
11251902	Leg Bag 600 ml, standard, 30 cm inlet, T-valve, soft backing	N	25	100
11251502	Leg Bag 600 ml, standard, 50 cm adjustable inlet, T-valve, soft backing	N	25	100
11421502	Leg Bag 800 ml, standard, 50 cm adjustable inlet, small lever , soft backing	N	25	100
11451502	Leg Bag 800 ml, standard, 50 cm adjustable inlet, T-valve, soft backing	N	25	100
11701101	Leg bag 1500 ml, 90 cm inlet tube, extra connctor, no outlet, soft backing	N	25	100
11751101	Leg bag 1500 ml, 90 cm inlet tube, extra connctor, free-flow, soft backing	N	10	100
12251902	Leg Bag 600 ml, standard, 30 cm inlet, fixed connector, sample port, T-valve, soft backing	S	10	100
12221502	Leg Bag 600 ml, standard, 50 cm inlet, fixed connector, small lever, soft backing	S	10	100
11401803	Leg Bag 750 ml, economy, 10 cm inlet, no outlet	N	–	100
11451803	Leg Bag 750 ml, economy, 10 cm inlet, T-valve	N	–	100
11451503	Leg Bag 750 ml, economy, 50 cm inlet, T-valve	N	–	100