



Product description

Order number	86.1180
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 104 x 10 mm, material: LD-PE, transparent, droplet size: 17 μ l - 22 μ l, 750 piece(s)/case

Product characteristics

Graduation	no
------------	----

Size

Sample volume	1 ml
Width of product	10 mm
Length of product	104 mm
Drop size	17 μ l - 22 μ l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Transfer pipette, 1 ml, (LxW): 104 x 10 mm, LD-PE,
transparent**

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

Packaging

Minimum order qty. 3000
Type of smallest case
subpackaging
Piece(s) / inner box 750
Piece (s) / outer case 3000
Piece(s) / pallet 252000
Depth of box 261 mm
Width of box 126 mm
Height of box 121 mm
Depth of case 285 mm
Width of case 260 mm
Height of case 260 mm
Case volume 0.0193 cbm
Weight of product 0.0005 kg
Weight of case 2.07 kg
EAN of inner box 4038917085649
EAN case 4038917004282

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-05-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT