



Product description

Order number 93.1650.020
Product description Insertion plate, Version 2, material: PP, 1 piece(s)/case

Product characteristics

Suitable for Safety Tray (93.1650) & e.g. 1x Multi-Safe vario 2000 (77.3893.020/.120), 1x rack (93.851 - 93.851.1xx), 2x dispenser boxes (95.813.001), 6x multi-purpose cup (75.1355), 1x disinfection spray

Size

Width of product 272 mm
Length of product 420 mm

Material & colours

Product material Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1440
Case volume	0.0009 cbm
Weight of product	0.162 kg
Weight of case	0.3 kg
EAN of inner box	4038917132701
EAN case	4038917076722

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

