



### **Product description**

|                     |                                                                                                                    |
|---------------------|--------------------------------------------------------------------------------------------------------------------|
| Order number        | 93.851.171                                                                                                         |
| Product description | Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, red, with handle, material: recycled PP |

### **Product characteristics**

|                        |        |
|------------------------|--------|
| Number of wells        | 50     |
| Storage capacity       | 50     |
| Format                 | 5 x 10 |
| Autoclavable (max. °C) | 121 °C |

### **Size**

|                     |        |
|---------------------|--------|
| Diameter of opening | 17 mm  |
| Width of product    | 109 mm |
| Height of product   | 45 mm  |
| Length of product   | 246 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10, red,  
with handle**

**Page 2**

**Material & colours**

|                   |                                      |
|-------------------|--------------------------------------|
| Product material  | recycled Polypropylene (recycled PP) |
| Colour of product | red                                  |

**Purity & certification**

|                  |               |
|------------------|---------------|
| Product category | not regulated |
| CE certified     | no CE mark    |
| Batched          | no            |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | -             |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 12            |
| Piece(s) / pallet             | 840           |
| Depth of case                 | 453 mm        |
| Width of case                 | 258 mm        |
| Height of case                | 165 mm        |
| Case volume                   | 0.0193 cbm    |
| Weight of product             | 0.1983 kg     |
| Weight of case                | 2.79 kg       |
| EAN of inner box              | 4038917121835 |
| EAN case                      | 4038917074971 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-05-09

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**