

## evercare Kraken Spike, Air- particle filter Blue with NFC

REF: 13493

Kraken<sup>®</sup> Spike for withdrawal and preparation of solutions, drugs and toxic substances. Minimizes risk for bacterial transfer, reduces risk for leakage and contamination of aerosols. Different version available for fulfilling all user needs. Kraken<sup>®</sup> Spike comes in three different colors with different specifications according to needs.

### Key Features

- ▶ Vented spike for withdrawing fluids from vials or semi-rigid containers.
- ▶ Easy handling due to ergonomic design.
- ▶ Low penetration force.
- ▶ Reduces contamination through safety snap cap.
- ▶ Multi-withdrawal possibility with a needle-free connector.
- ▶ Available with an air filter or a combination of air and fluid filters.
- ▶ Red Kraken Spike for cytostatic drugs.
- ▶ PVC-free, DEHP-free, Latex-Free



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
13493	evercare Kraken Spike, Air- particle filter blue with NFC, Blue	0,1µm venting & 5µm particle filter	100 / 1000

Sales Unit: Transport Carton

## Regulations

N/A



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-01  
EN ISO 10993-7

**Compliance with legal requirements:** Medical Devices Directive 93/42/EEC, MDD Class Is

**Compliance with product standards:** EN ISO 80369-7

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 20417, EN ISO 15223-1

## Material

Product Part	Raw Material & Characteristics
<b>Green &amp; Blue Spike Red Chemo Spike Spike Shield Spike cover Air Filter, Green, Blue &amp; Red Spike Fluid Filter, Blue &amp; Red Spike Green, Blue &amp; Red Spike Needle Free Connector Snap Cap</b>	Acrylonitrile Butadiene Styrene (ABS) Copolymer ETT Acrylonitrile Butadiene Styrene & Copolymer Polypropylene (PP) PTFE Membrane PET Colorant Polypropylene Cosmoplene & Silicon Polypropylene (PP)

## Warnings & Recommendations

### Warnings

**Disposal Recommendation** Dispose in accordance with local regulations.

**Shelf Life** 5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
Single Pack	Medical grade paper, PP/PE film	105 x 55
Dispenser Pack	Cardboard	215 x 160 x 205
Transport Carton	Corrugated cardboard	280 x 225 x 405
Packing Tape	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
13493	8902095134931	38902095134932	58902095134936

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.