

Surgical Cap Lady

REF: 2794-01

Surgical Cap Lady is a surgical cap that has been specially developed to provide good coverage of the hair. Size can be adjusted via a binding function in the neck. The surgical cap is made in a soft and breathable nonwoven material.

The use of headwear is proven to reduce the risk of particles falling into the surgical wound.

Key Features:

- ▶ Single use surgical cap for coverage of long hair
- ▶ Soft viscose material, top of cap made of thinner nonwoven material for better breathability
- ▶ Tied in back
- ▶ Color: green/white
- ▶ Universal size



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2794-01	Surgical cap Lady, green	Universal size	100/600

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices

[**NATURAL LATEX FREE**]    

REACH regulation: The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-10

Compliance with legal requirements: Medical Devices Directive 93/42/EEC, Class I

Compliance with product standards: N/A

Quality standards: EN ISO13485

Label information and packaging compliance: EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristic
--------------	-------------------------------

Main material	Viscose nonwoven, green
---------------	-------------------------

Top	Spunbond PP nonwoven, white
-----	-----------------------------

Sewing thread	Nylon, white
---------------	--------------

Warnings and Recommendations

Warnings: The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

Storage Recommendation Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.

Disposal Recommendation Product and packaging materials to energy recovery or mixed waste.

Shelf Life 3 Years

Packing Information

Package	Material	Size (L x W x H) mm
---------	----------	-----------------------

Single Pack	N/A	N/A
-------------	-----	-----

Dispenser Pack	Carton, TCF quality	300 x 185 x 145 mm
Transport Carton	Carton, TCF quality	385 x 310 x 455 mm
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2794-01	N/A	 6 438129 319811	 6 438129 519815

ONEMED

Manufactured in China for
OneMed Group Oy
 Address: Metsäläntie 20
 Postal Address: P.O. Box 10
 FI-00321 Helsinki, Finland
 Phone +358 20 786 6800
www.onemed.com
 Domicile Helsinki
 Business ID 2039640-1