



Short-style pediatric mask 001261



Under-the-chin-style pediatric mask 001263



Adult high-concentration mask 001207

AirLife™ Aerosol Masks

These soft, high-quality AirLife™ disposable aerosol masks are made from high-grade resins. The smooth, feathered edges are designed for patient comfort while reducing irritation points. All AirLife™ masks are designed for anatomical compatibility. A full line is available to meet requirements in pediatric and adult sizes.

Cat. No.	Description	Qty.
Pediatric Masks		
001261	Short-style mask	50/cs
001263	Under-the-chin-style mask	50/cs
Adult Masks		
001206	Mask only, under-the-chin style	50/cs
001207	Mask only, high concentration	50/cs
001209	Mask with 5' (60") corrugated tubing	50/cs
001208	Mask with 6' (72") corrugated tubing	50/cs



Tracheostomy mask 001225

AirLife™ Tracheostomy Masks

The AirLife™ disposable tracheostomy masks are available in pediatric and adult sizes, and incorporate several patient features. The neckband is made from a comfortable, nonbiting material, and the swivel tubing connector permits access from either side of the patient. The special strap clips allow the mask to be removed with minimal disturbance to the patient.

Cat. No.	Description	Qty.
001226	Pediatric mask	50/cs
001225	Adult mask	50/cs
001227	Adult mask, large	50/cs



Face tent mask 001220

AirLife™ Face Tent Masks

The AirLife™ disposable face tent masks are available in an adult under-the-chin style to be used with large bore aerosol tubing for high humidity aerosol therapy. The face tent is molded of soft vinyl. The under-the-chin design and clarity offer patient comfort, clear vision and ease of speech. The face tent incorporates a thick wall to give more stability to the mask and is available with 5' (60") or 6' (72") corrugated tubing for greater versatility.

Cat. No.	Description
001220	Mask only
001222	With 5' (60") corrugated tubing
001221	With 6' (72") corrugated tubing