



Product description

Order number 85,116
Product description Needle for S-Monovette® 20G x 1 1/2

Product characteristics

Colour code Yellow
Design Needle
Safety product no
Connection Membrane screw cap

Size

Diameter of the needle (external) 20G/0,9 mm
Length of the product 38 mm
Length of the needle 38 mm

Purity & certification

Product in accordance with MDD or IVDD Medical Product
CE certified CE0197
Sterilisation Gamma irradiation
Pyrogen-free Pyrogen-free
Lot or serial number always batched

Packaging

Type of packaging individual sterile packaging
Pcs. / innerbox 100
Pcs. / box 1000
Pcs. / pallet 78000
Minimum order qty.(pcs.) 1000
Depth of case 368
Width of the case 233
Height of the case 290
Case volume 0,0249 cbm
Weight of the case 3,02 KG
EAN of the inner box 4038917085519
EAN of the outer box 4038917004169

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

