



Product description

Order number	86,1175
Product description	Transfer pipette 6 ml, without graduation, Length: 146 mm, Drop size: 20 - 30 µl

Product characteristics

Graduation	No
------------	----

Size

Sample volume	6 ml
Max. volume	6 ml
Length of the product	146 mm
Drop size	20-30 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	Transparent

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	not always batched

Packaging

Pcs. / innerbox	400
Pcs. / box	1600
Pcs. / pallet	44800
Minimum order qty.(pcs.)	1600
Depth of case	533
Width of the case	333
Height of the case	250
Case volume	0,0444 cbm
Weight of the case	3,31 KG
EAN of the inner box	4038917085618
EAN of the outer box	4038917004251

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

