



## Product description

Order number	72.690.475
Product description	Micro tube 1.5 ml with attached lid with a wide hinge, PP, blue, 500 pieces packed in the bag

## Product characteristics

Graduation	Yes
Cap	Attached
Recommended centrifugation	20000 x g

## Size

Sample volume	1,5 ml
Diameter	10,8 mm
Length of the product	39 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	Blue
Cap material	Polypropylene (PP)
Colour of cap	Blue

## Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	not always batched

## Packaging

Pcs. / innerbox	500
Pcs. / box	5000
Pcs. / pallet	180000
Minimum order qty.(pcs.)	5000
Depth of case	398
Width of the case	398
Height of the case	300
Case volume	0,0475 cbm
Weight of the case	5,69 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification  
Micro tube 1.5ml, blue**

**Page 2**

EAN of the inner box            4038917161688  
EAN of the outer box            4038917015226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2015-11-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**