



### **Product description**

Order number 20.1288.100  
Product description Microvette® 200 µl, K3 EDTA, violet US code, print

### **Product characteristics**

Type of collection Capillary  
Type of preparation K3 EDTA  
Colour code US  
Label/ Print With printed label  
Colour of print/label Violet  
Cap Quadruple-threaded cap  
Application Capillary blood collection  
Base shape Flat base  
Recommended centrifugation 2000 x g  
Centrifugation temperature Room temperature  
Centrifugation time 5 min  
Shape of the inner tube Cylindrical

### **Size**

Sample volume 200 µl  
Diameter 10,8 mm  
Length including cap 47,6 mm  
Length excluding cap 46,6 mm

### **Material & colours**

Product material Polypropylene (PP), Polyvinylchloride (PVC)  
Colour of product Transparent  
Cap material High Density Polyethylene (HD-PE)  
Colour of cap Violet

### **Purity & certification**

Product in accordance with In-Vitro Product  
MDD or IVDD  
CE certified CE  
Lot or serial number always batched

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





## **Packaging**

Pcs. / innerbox	100
Pcs. / box	1000
Pcs. / pallet	121000
Minimum order qty.(pcs.)	1000
Depth of case	338
Width of the case	228
Height of the case	160
Case volume	0,0123 cbm
Weight of the case	3,28 kg
EAN of the inner box	4038917127523
EAN of the outer box	4038917070393

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

