



Product description

Order number	75.1337.500
Product description	Medicine cup 30 ml, PP, natural, 75 pieces packed in the PE bag

Product characteristics

Graduation	Yes
Example of use	Collection tube for urine samples

Size

Sample volume	30 ml
Diameter	37 mm
Height of the product	40 mm
Length of the product	41 mm
Length excluding cap	40 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent

Purity & certification

Product in accordance with MDD or IVDD	Medical Product
CE certified	CE
Lot or serial number	not always batched

Packaging

Pcs. / innerbox	3750
Pcs. / box	3750
Pcs. / pallet	180000
Minimum order qty.(pcs.)	3750
Case volume	0,0319 cbm
Weight of the case	6 KG
EAN of the inner box	4038917127721
EAN of the outer box	4038917043144

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

