

# Connecting Tubes for Drainage Bags

We offer a wide range of connecting tubes for all types of drainage bags. The connecting tubes are available with or without stopcocks.

- Flexible, kink resistant tubing.
- Clinical versatility.

## Connecting Tubes, Adapters and Stopcocks Box of 25

- 4039 43 015** 15 cm female-olive
- 4040 17 015** 15 cm male-funnel
- 4040 17 030** 30 cm male-funnel
- 4040 18 015** 15 cm one-way stopcock/male-funnel
- 4040 18 030** 30 cm one-way stopcock/male-funnel
- 4040 22 000** Male/three-way stopcock/funnel
- 4040 24 000** Male/funnel
- 4040 40 015** 15 cm male - male luer lock
- 4040 40 030** 30 cm male - male luer lock
- 4040 43 015** 15 cm male -3-way stopcock/male luer lock
- 4040 43 030** 30 cm male -3-way stopcock/male luer lock
- 4043 17 050** 50 cm male - tube
- 4043 43 015** 15 cm male luer lock/3-way stopcock - funnel
- 4043 43 030** 30 cm male luer lock/3-way stopcock - funnel
- 4047 43 015** 15 cm male luer lock - olive
- 8011 00 102** Stopcock – low pressure, max pressure 50 psi
- 8011 00 103** 3-way stopcock – low pressure, max pressure 50 psi



4040 17



4040 22



4040 24



4040 40




4040 43



4043 17



4047 43



8020 20

**8020 20 000** Female/olive with stopper

8020 21

**8020 21 000** Male/olive with stopper

Box of 10



**Angiotech Denmark**

Knud Bro Alle 4 • DK-3660 Stenløse

customer@angio.com • www.AngioEduPro.com