



## evercare® Surgical cap Tina

REF: 2782-01, 2783-01

evercare® Headweares are designed to help prevent foreign particles such as hair or skin scales falling from medical personnel into the wound during surgical interventions. Our surgical caps are suitable for most interventions and situations and are available in a wide range of colors and designs. evercare® surgical cap Tina serves staff who wants maximum hair coverage. Made in a soft and breathable viscose nonwoven material (100% fossil-free). There is a good coverage also for long hair and an elastic in the back for a secure fit. Available in two colors.

### Key Features

- ▶ Reduces the risk of particles falling into the wound
- ▶ Universal size
- ▶ Designed for maximum hair coverage
- ▶ Elastic in the back




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2782-01	Surgical cap Tina, Green	Universal	100 / 600
2783-01	Surgical cap Tina, Blue	Universal	100 / 600

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
Main material	Viscose nonwoven
Elastic band	Polyester yarn and Synthetic rubber, latex free, white
Sewing thread	Nylon, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	300 x 185 x 145
<b>Transport Carton</b>	Carton, TCF or ECF quality	385 x 310 x 455
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>2782-01</b>	N/A	6438129319750	6438129519754
<b>2783-01</b>	N/A	6438129319767	6438129519761

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.