



evercare® Surgical cap Tristan

REF: 2791-01

evercare® Headweares are designed to help prevent foreign particles such as hair or skin scales falling from medical personnel into the wound during surgical interventions. Our surgical caps are suitable for most interventions and situations and are available in a wide range of colors and designs. evercare® surgical cap Tristan is a popular design in the USA. Its smaller design does not fully cover ears and long hair. Made in a soft and breathable viscose nonwoven material (100% fossil-free). The top of the cap is made of thinner nonwoven material (white) for better breathability. The tie band in the back allows for individual fit.

Key Features

- ▶ Reduces the risk of particles falling into the wound
- ▶ Universal size
- ▶ Top of cap in thinner nonwoven material for better breathability
- ▶ Tie band in the back allows for secure fit




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2791-01	Surgical cap Tristan, Green	Universal	100 / 600

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Body	Viscose nonwoven, green
Cap	Spunbond PP nonwoven, white
Tie bands	Spunbond PP nonwoven, white
Sewing thread	Nylon and polyester, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	300 x 185 x 145
Transport Carton	Carton, TCF or ECF quality	385 x 310 x 455
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2791-01	N/A	6438129319804	6438129519808

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.