

PRODUKTINFORMATION

GRIPPER PLUS[®] POWER P.A.C



KANYL TILL SUBCUTAN VENPORT

- Klarar 300psi eller max 5ml/sek.
- Ett avtagbart grepp säkerställer en kontrollerad nålplacering.
- En icke absorberande skumdyna för ökad patientkomfort.
- Integrerad anslutnings slang ca 20cm med luerkoppling och färgkodad avstängningsklämma.
- Nålskydd för att minska risken för stickskador.
- Ger en bra patientkomfort samtidigt som den är lätt och säker att avlägsna.
- Finns med eller utan nålfri Y-koppling, NAC.



Benämning	Art. nr	Art nr. NAC	Storlek	Längd	Enhet	Avdfp
Gripper Plus [®] Power P.A.C.	21-3368-24	21-3468-24	19G	19mm	st	12
Gripper Plus [®] Power P.A.C.	21-3364-24	21-3469-24	19G	25mm	st	12
Gripper Plus [®] Power P.A.C.	21-3365-24	21-3470-24	19G	32mm	st	12
Gripper Plus [®] Power P.A.C.	21-3369-24	21-3464-24	20G	16mm	st	12
Gripper Plus [®] Power P.A.C.	21-3367-24	21-3465-24	20G	19mm	st	12
Gripper Plus [®] Power P.A.C.	21-3362-24	21-3466-24	20G	25mm	st	12
Gripper Plus [®] Power P.A.C.	21-3363-24	21-3467-24	20G	32mm	st	12

Smiths Medical Sverige AB

Box 1143, 164 22 Kista
 Tlf. 08 - 594 772 50
 Fax.08 - 594 772 69
 info.sweden@smiths-medical.com
www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol[®] indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
 ©2012 Smiths Medical. All rights reserved. SE_14gGPI_GPP-1