

Rack for Ø 30 mm tubes



Technical modifications reserved

- 18 Apertures for 50 ml centrifuge tubes or 25 ml S-Monovettes of the Formalin System
- Autoclavable
- Convenient handles
- Dimensions in mm (W x D x H): 295 x 114 x 66

Ordering Information

Order No.	Description	Packaging
93.1707	Rack, blue, for Ø 30 mm tubes, 18 apertures	1/bag 10/case

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-122
info@sarstedt.com
www.sarstedt.com

 **SARSTEDT**