


## High Performance Breathing Circuits for Optiflow™



The RT-series breathing circuits were designed specifically for use with Fisher & Paykel Healthcare humidification systems. Spiral heater wire technology enables superior condensate reduction performance in varying environments by consistent distribution. Heater wire performance is uniquely keyed to humidifier algorithms. Dual-spiral heater wire ensures even heat distribution leading to reduced condensate.

### Adult Optiflow Breathing Circuits

Part No	Description	Quantity
RT202	Single Limb Adult Breathing Circuits Kit	10/box
RT102	Single Limb Adult Breathing Circuit (excl. Chamber)	10/box
RT241	Single Limb Adult Breathing Circuit Kit for use with MR880 humidifier	10/box
900PT501	Patient Circuit and Chamber Kit for use with AIRVO/myAIRVO humidifier	10/box