

Optiflow Junior Oxygen Tubing New Release

Fisher & Paykel Healthcare is pleased to announce the release of the Optiflow Junior Oxygen Tubing.



This oxygen tubing is designed to connect the Optiflow Junior nasal cannula to a portable oxygen source. This allows the delivery of oxygen as an interim measure during transport, additionally it can be used as a standard oxygen tubing.

Flows up to 2 L/min can be provided to patients. This tubing has been designed to be used on a single-use basis.

Each Optiflow Junior Oxygen Tube is individually packed and distributed as a box of ten.

Relevant Part Number

These can be ordered using the OPT014 part number.

Please contact your local representative for further information about this or any other Fisher & Paykel Healthcare product or therapy, or visit www.fphcare.com.