

Optiflow Junior Adaptors Kit

New Release

Fisher & Paykel Healthcare is pleased to announce the release of the Optiflow Junior Adaptors Kit. This kit is designed to enable aerosolized medications to be delivered through the Optiflow Junior cannulae and gas sampling using the Luer lock. Insertion of the nebulization adaptor between the patient end of the circuit and the interface enables the addition of aerosolized medications to the flow of gas.

Each of the Optiflow Junior Adaptor Kits is for single use and contains:

- An in-line nebulization adaptor with both a 22mm male port and a luer port.
- A 22mm elbow, female to male.

The Optiflow Junior Adaptors Kit contains a pair of the two connectors which are packed and distributed as a box of ten pairs.

Relevant Part Number

These kits can be ordered using the OPT016 part number.



Please contact your local representative for further information about this or any other Fisher & Paykel Healthcare product or therapy, or visit www.fphcare.com.