



Alaris™ VP Series DEHP FREE infusion set with 15µm filter 70593

- Set Length 300 cm
- Priming Volume 27 ml
- Bore 3 mm
- 1 Spiked Drip Chamber(s) (Vented/Unvented)
- Drops per ml 20
- Hydrophobic Air Inlet Filter (with Cap)
- 15 µm Drip Chamber Filter
- 1 Roller Clamp(s)
- Alaris™ Safety Clamp
- Silicone pump segment
- Pressure Sensor Disc
- 1 Male Luer(s)
- Polyvinyl Chloride (PVC)
- Natural rubber latex is not part of the material formulation
- DEHP is not part of the material formulation
- Quantity 80

