

# droplet<sup>®</sup>

personal lancets



CE 0344

droplet<sup>®</sup> personal lancets

ultra thin and universal

# droplet®

personal lancets

Comfort  
MAX

33G

- ComfortMax – new thinnest needle and highest comfort
- delicate – silicone-coated needle
- safe – gamma-sterilized
- universal – fits most lancing devices



**Droplet® lancing device - specially designed for highest comfort in your daily blood testing.**

- personalized – 5 lancing penetration depths
- safe – lancing ejection system to allow safe removal of used lancets
- universal – compatible with most popular lancets
- gentle – smooth puncture with no vibrations, now even greater comfort with Droplet® lancets 33G

**Optimal use with Droplet® lancets**



30G


ultra  
thin

- ultra thin and precise
- delicate – silicon-coated needle minimizes pain
- especially recommended for patients with delicate skin, including children
- universal – fits most lancing devices available on the market



**Results of multiple use:**

- high risk of infection
- increase of pain
- tissue damage
- slow healing



28G

- thin and perfectly ground
- recommended for adults and patients with thicker skin
- universal – fits most lancing devices available on the market



1x use

**Remember!**

Use Droplet® lancets **only once** and dispose of after use

