

# C-A-T® Combat Application Tourniquet

Gen 7, EMS

---



---

Item No:	CL-0002
Manufacturer:	C-A-T Resources LLC
Certificates:	CE
Classification:	Class I
Size:	95.25cm / 37.5Inch
Color:	Orange
Country of origin:	US
NSN no:	6515-01-612-7129
Tariff:	90189084
Weight:	85g
Dimensions:	17 x 4 x 4cm
Shelf life:	N/A

A Combat Application Tourniquet makes a significant difference in a life or death situation as it is designed to stop bleedings on extremities.

The C-A-T tourniquet design offers the same benefits to EMS users as it has over the past years for America's military soldiers. With its Velcro and safety buckle, and correct training the C-A-T Tourniquet can easily be operated by one hand only.

The C-A-T Tourniquet is made of a special composite material. Latex free. C-A-T is approved by CoTCCC.