



# Reusable cuffs that stand the test of time

Philips blood pressure solutions are simple to use, perform as expected, and save time. These reusable cuffs are easy to care for with a broad line of cleaning agents, helping to prevent the spread of hospital-acquired infections.

## Easy Care for an easier day

With rounded edges for greater comfort, and designed for use in most areas of the hospital, Easy Care Multi-Patient Cuffs are durable, come in a wide range of sizes, and are treated with Micropel™, an antimicrobial agent that prevents the growth of fungal, bacterial, and algal organisms on the cuff.<sup>1</sup>



## Easy Care Multi-Patient Cuffs

Single-hose cuff	Double-hose cuff	Cuff size	Color
M4552B	M4562B	Infant	Orange
M4552B5			
M4553B	M4563B	Pediatric	Green
M4553B5			
M4554B	M4564B	Small adult	Light blue
M4554B5			
M4555B	M4565B	Adult	Navy
M4555B5			
M4556B	M4566B	Adult extra-long*	Navy
M4556B5			
M4557B	M4567B	Large adult	Burgundy
M4557B5			
M4558B	M4568B	Large adult extra-long**	Burgundy
M4558B5			
M4559B	M4569B	Thigh	Grey
M4559B5			

Accuracy of extra-long cuffs is dependent on proper application according to cuff labeling.

\*26.0 inches (66.0 cm) total cuff length.

\*\*32.9 inches (83.5 cm) total cuff length.

<sup>1</sup>The mechanism of action of this agent does not contribute to the development of drug-resistant microbial strains. Reduction in colonization or microbial growth on the device has not been shown to correlate with a reduction in infections in patients. Clinical studies to evaluate reduction in infection have not been performed.

**PHILIPS**  
sense and simplicity