

# Counting chambers, C-Chip

Leverantör: NANOENTEK









Quartz grade optical plastic counting chambers applicable in blood analysis (hematology).

- ▶Guaranteed repeatability and reliability, eliminate the need for washing and reuse
- ▶Disposable hemocytometers that do not need coverslips
- ▶Chip to reduce exposure to infectious samples and contamination

Suitable for cell counting and culture, concentration measurement and viability test, in bacteria and fungal spore counting, IVF, IUI tasks of sperm counting.

## ORDER

VWR art.nr	Antal	Benämning	Typ	Depth of chamber	B×D	Förp.
631-1094		Counting chambers	Bürker	0,1 mm	25×75 mm	50
631-1095		Counting chambers	Bürker-Türk	0,1 mm	25×75 mm	50
631-1098		Counting chambers	Imp. Neubauer	0,1 mm	25×75 mm	50
NANODHC-N04		Counting chambers	Neubauer	0,1 mm	25×75 mm	50
631-1096		Counting chambers	Fuchs-Rosenthal	0,2 mm	25×75 mm	50
631-1097		Counting chambers	Malassez	0,2 mm	25×75 mm	50