Pulmodyne®



BiTrac MaxShield™

Fixed and Select Options to Suit a Variety of Patient Needs

- MaxShield comes in a wide variety of sizes and elbow options.
- Assists in management of claustrophobia.
- MaxShield comes in a Fixed line and a Select line with interchangeable elbows.
- A quick and easy fit for patients.

BiTrac MaxShield Select™ with Interchangeable Anti-Asphyxia Elbow

Part No.	Sizes	Case/Qty.
313-9559	Adult Extra Large	10
313-9558	Adult Large	10
313-9557	Adult Small	10

BiTrac MaxShield Select™ with Interchangeable Anti-Asphyxia with Leak Elbow

Part No.	Sizes	Case/Qty.
313-9555	Adult Extra Large	10
313-9552	Adult Large	10
313-9554	Adult Small	10

BiTrac MaxShield SE Select™ with Interchangeable Standard Elbow

Part No.	Sizes	Case/Qty.
313-9556	Adult Extra Large	10
313-9553	Adult Large	10
313-9551	Adult Small	10

Interchangeable Elbows

Part No.	Sizes	Case/Qty.
313-9599	Standard Elbow (dual limb circuit)	10
313-9593	Anti-Asphyxia with Leak Elbow (single limb circuit)	10
313-9597	Standard Bronchoscopy Elbow (dual limb circuit)	10

BiTrac MaxShield™ with fixed Anti-Asphyxia Elbow

Part No.	Sizes	Case/Qty.
313-9057W	Adult Extra Large	10
313-9052W	Adult Large	10
313-9054W	Adult Small	10
DiTues CE MayCh	delating one recording	

BiTrac SE MaxShield™ with fixed Standard Elbow

Part No.	Sizes	Case/Qty.
313-9056W	Adult Extra Large	10
313-9053W	Adult Large	10
313-9055W	Adult Small	10















Pulmodyne®

For more information about the BiTrac MaxShield Select™, BiTrac MaxShield™, and our complete line of emergency airway products visit our website at

2055 Executive D Indianapolis, IN 46241



....bringing change to life.®

pulmodyne



